

Cause-and-effect analysis as a tool to improve the reproducibility of nanobioassays: four case studies

Elijah J. Petersen*¹, Cordula Hirsch², John T. Elliott¹, Harald F. Krug³, Leonie Aengenheister², Ali Talib Arif^{4, 5}, Alessia Bogni⁶, Agnieszka Kinsner-Ovaskainen⁶, Sarah May², Tobias Walser⁷, Peter Wick², Matthias Roesslein²

1) Biosystems and Biomaterials Division, Material Measurement Laboratory, National Institute of Standards and Technology (NIST), Gaithersburg, MD 20899

2) Empa, Swiss Federal Laboratories for Material Testing and Research, Particles-Biology Interactions Laboratory, CH-9014 St. Gallen, Switzerland

3) NanoCASE GmbH, St. Gallerstr. 58, 9032 Engelburg, Switzerland

4) Institute for Infection Prevention and Hospital Epidemiology, University Medical Center Freiburg, University of Freiburg, Faculty of Medicine, University of Freiburg, D- 79106 Freiburg, Germany

5) Kurdistan Institution for Strategic Studies and Scientific Research (KISSR), Qirga - Sulaimani, Iraq


6) European Commission, DG Joint Research Centre, 21027 Ispra, VA, Italy

7) Vereala GmbH, Eyhof 34, 8047 Zürich, Switzerland

*Corresponding author: **E.J.P.**, elijah.petersen@nist.gov, Tel: (301)-975-8142

19 **TOC graphic**

Sources of variability in assays for use with nanoparticles



20

21 **Abstract**

22 One of the challenges in using *in vitro* data to understand the potential risks of engineered nanomaterials
23 (ENMs) is that results often differ or are even contradictory among studies. While it is recognized that
24 numerous factors can influence results produced by nanobioassays, there has not yet been a consistently
25 used conceptual framework to identify key sources of variability in these assays. In this paper, we use
26 cause-and-effect analysis to systematically describe sources of variability in four key *in vitro*
27 nanobioassays: the DCF (2',7'-dichlorofluorescein) assay, an enzyme-linked immunosorbent assay (ELISA)
28 for measuring interleukin-8, a flow cytometry assay (Annexin V/Propidium Iodide), and the Comet assay.
29 These assays measure endpoints that can occur in cells impacted by ENMs through oxidative stress, a
30 principle mechanism for ENM toxicity. The results from this analysis identify control measurements to test
31 for potential artifacts or biases that could occur during conduct of these assays with ENMs. Cause-and-
32 effect analysis also reveal additional measurements that could be performed either in preliminary
33 experiments or each time the assay is run to increase confidence in the assay results and their
34 reproducibility within and among laboratories. The approach applied here with these four assays can be
35 used the support the development of a broad range of nanobioassays.

36

37

38 Introduction

39 Engineered nanomaterials (ENMs) often have unique or enhanced properties such as high surface
40 reactivity and quantum confinement compared to bulk materials of the same elemental composition.
41 These properties can be utilized in a broad range of commercial applications such as environmental
42 remediation, biomedical application, textiles, and renewable energy.¹⁻⁵ However, ENMs may also pose
43 potential risks to human health or the environment as a result of these same properties since a fraction
44 of the ENMs will be intentionally or unintentionally released during the production and life cycle of these
45 products.⁶⁻¹⁵ Many assays used to assess biological effects of dissolved chemicals may require
46 modifications prior to use with ENMs because these materials often behave differently in the test media
47 (e.g., agglomeration or dissolution) and may cause artifacts in the test results.¹⁶⁻²²

48 The most frequently used tests to assess the potential hazards of ENMs are *in vitro* toxicity
49 assays.¹⁷ They have the advantage of enabling higher throughput testing which is not possible when
50 conducting tests with larger and more complex organisms (e.g., rats or fish). High throughput testing can
51 be used to develop categories for ENMs based on results from these assays, the ENM composition, and
52 physicochemical properties.²³⁻²⁶ In some cases, there is evidence that *in vitro* assay results might also be
53 used to estimate the toxicity observed with more complex organisms such as pulmonary toxicity with
54 rats.^{27,28} In addition, *in vitro* tests such as cytotoxicity assays can be used in a tiered testing system to
55 assess potential hazards caused by ENMs. This can support the prioritization of additional testing using
56 larger organisms thus reducing the number of animals used.²⁷⁻²⁹

57 One of the most substantial limitations for *in vitro* testing with ENMs is a lack of robust,
58 “validated” methods for these materials.^{27,30-33} The challenges of assessing ENMs with existing test
59 methods has also been recognized by the test guideline program of the Organization for Economic Co-
60 operation and Development (OECD). The OECD launched the testing program of ENMs in 2007 to ensure
61 that the approaches for hazard, exposure, and risk assessment for ENMs are science based and
62 internationally harmonized.³⁴ One objective of this program is to explore the potential application of
63 alternative methods for testing of ENMs. Validated methods have been developed for dissolved chemical
64 substances, yet additional work is needed in some cases to adapt these methods for use with ENMs given
65 their different behaviors in the assay system. For example, there is a possibility for ENMs to cause artifacts
66 in these assays such as by adsorbing reagents of the test system,³⁵⁻³⁸ damaging biomolecules after a cell
67 assay has concluded (e.g., ENM induced photoactivity as a result of laboratory lighting),^{21,39} or producing
68 a signal (e.g., fluorescence or absorbance) similar to the measurand for the assay.^{19,20} Biases caused by
69 these artifacts and other important, yet often uncontrolled, aspects of the assay such as cell-seeding
70 density may be the source of contradictory results on the toxicity of various ENMs. Another complication
71 for the development of robust assays for ENMs is the dispersion protocol, which, if not standardized, may
72 lead to different ENM suspensions even when using the same starting material.^{17,18,40-42} Testing the ENMs
73 with same composition but from different manufacturers may generate different results, hindering
74 comparison among *in vitro* assay results even within a single laboratory.⁴³

75 One approach that we have used previously to comprehensively evaluate potential sources of
76 uncertainty for the 3-(4,5-dimethylthiazol-2-yl)-5-(3-carboxymethoxyphenyl)-2-(4-sulfophenyl)-2H-

77 tetrazolium (MTS) cell viability assay is cause-and-effect (C&E) analysis.¹⁵ C&E analysis is an approach,
78 originally utilized in quality manufacturing and propagation of measurement error in analytical
79 chemistry.⁴⁴ It highlights through C&E diagrams, prepared based on expert judgment, aspects of the assay
80 anticipated to most strongly influence the variability of the assay result. C&E diagrams, in particular,
81 provide a graphical representation of potential sources of variability in assays. These diagrams can then
82 be used to add process control measurements to an assay protocol to track key potential sources of
83 variability.¹⁵ Understanding these sources of variability and the robustness of the assay to minor,
84 unintended changes in the protocol (e.g., temperature, stability of reagents, time range for incubation)
85 enables scientifically informed, instead of *ad hoc*, choices about components of the assay to monitor, and
86 allowable ranges for different steps of the assay (e.g., temperature or incubation duration). This supports
87 reducing interlaboratory variability during interlaboratory testing and eventual approval for usage in
88 regulatory testing.⁴⁵ For example, during sensitivity testing of the MTS assay, we determined that the cell
89 pipetting caused the highest variability in the assay result.¹⁵ These findings and the use of process control
90 measurements supported successful interlaboratory testing of the protocol where outlier results from
91 one laboratory could be traced to a different interpretation of a single assay step and corrected by revising
92 the procedure.⁴⁶ However, it was unclear from this case study about how to apply this approach to other
93 nanocytotoxicity assays and if the same control measurements would be sufficient to capture the most
94 important steps in these assays, particularly if the assays required more complex steps than the relatively
95 straightforward MTS assay.

96 To address this topic, a workshop took place in St. Gallen, Switzerland, during June 2015 entitled
97 “Cause-and-effect analysis: A new approach for developing robust nanobioassays” with 17 participants
98 from 4 countries. During the workshop, the application of C&E analysis was evaluated on four additional
99 cell-based assays, selected by the workshop participants, for use with ENMs: the DCF (2',7'-
100 dichlorofluorescein) assay, an enzyme-linked immunosorbent assay (ELISA) for measuring interleukin-8
101 (IL-8), a flow cytometry assay (Annexin V/Propidium Iodide (PI)), and the Comet assay. These assays can
102 be used to measure a possible cascade of events that can occur in cells impacted by ENMs through
103 oxidative stress, a condition in which the amount of intracellular ROS produced overwhelms the cells'
104 antioxidant defense capacities.^{47,48} Oxidative stress has been shown in numerous nanotoxicity studies to
105 be the principal mechanism causing toxicity given the potential for a broad range of ENMs to produce
106 reactive oxygen species (ROS) or activate redox reactions. Because of this, reactivity is among the most
107 promising intrinsic ENM properties used in ENM hazard categorization.^{47,49-52} Given that it is not yet clear
108 which combination of assays measuring oxidative stress is most predictive to estimate *in vivo* effects,⁴⁷
109 the previously mentioned four assays were evaluated: the DCF assay which directly measures intracellular
110 ROS, the ELISA which measures the increase in cytokines indicative of inflammation, the flow cytometry
111 assay which measures the mechanism of cell death which can be caused if the cells are exposed to
112 sustained oxidative stress, and the Comet assay which measures DNA damage, another potential outcome
113 of increased intracellular ROS levels. However, it should be noted that positive responses with these
114 assays can also be caused by mechanisms other than oxidative stress. For each assay, we performed a
115 C&E analysis and designed control measurements to assess for ENM interference with the assay. The
116 extent to which data from positive chemical control samples can provide insights into how well the assay
117 is functioning is also discussed. Lastly, next steps are described that could be taken with these assays to
118 more fully evaluate different expected sources of variability.

119

120 **Overview of commonly used in vitro-based cytotoxicity assays**

121 An overview for each of the four assays investigated during the workshop is provided in Table 1
122 and each assay will be briefly described. The full protocols for each assay are provided in the Supporting
123 Information (SI). Flow charts which describe each step of the assay are shown in Figure 1. More detailed
124 descriptions and protocols for each of these assays are available in the SI. The DCF assay is designed to
125 measure production of reactive oxygen species (ROS) by ENMs, one of the most common mechanisms of
126 ENM toxicity. The assay works by assessing the rate at which a non-fluorescent dye is chemically reduced
127 by ROS to form a fluorescent dye. This assay can be conducted under acellular or cellular conditions. In
128 the acellular version of the assay, there are additional steps to remove the diacetate moiety while
129 deacetylation occurs intracellularly in the cellular version of the assay. In the ELISA, the release of
130 cytokines, molecules that are indicative of inflammatory reactions, is quantified. The production of a
131 specific cytokine, in this example interleukin-8 (IL-8), is measured using the sandwich ELISA method which
132 works as follows: 1) adsorbing a primary antibody onto the surface of a high-affinity binding microwell
133 plate, 2) having the antibody bind the protein of interest in the cell culture supernatant, 3) adding a second
134 antibody to bind to the same protein of interest but at a different epitope, 4) adding horseradish
135 peroxidase linked to avidin to initiate an enzymatic reaction with tetramethylbenzidine (TMB) as the
136 substrate, and 5) quantifying the TMB using absorbance on a plate reader. In the third assay, flow
137 cytometry can be used to evaluate the quantity of apoptotic and necrotic cells. Cells are treated with
138 markers to distinguish between (i) viable cells, (ii) cells which are not viable but have the membrane intact
139 (apoptotic cells), and (iii) cells which are not viable and also have undergone membrane disintegration
140 (necrotic cells). The cells are subsequently analyzed using flow cytometry. In contrast to the *in vivo*
141 situation, apoptotic cells cannot be removed by tissue macrophages *in vitro*. Their membranes start to
142 disintegrate and cells stain double-positive for both markers (late apoptotic cells). Lastly, the purpose of
143 the Comet assay is to assess the potential genotoxicity of a compound through measuring the degree of
144 DNA damage that has occurred. The DNA integrity of cells is evaluated after gel electrophoresis, and the
145 length and quantity of DNA in the Comet tail relates to the extent of DNA damage.

146

147 **Control charting data**

148 Control charting data of the positive and negative control results (i.e., results for these control
149 measurements from assays performed on multiple different days) from the four assays obtained at Empa
150 using A549 cells are provided in Figures S1 through S4. The control charting data was also evaluated to
151 assess if there was a statistically significant correlation among the control charting data for different
152 components of data in the control charts (e.g., the negative control data and the positive control data)
153 (e.g., Figures S5 and S6). The nonparametric Spearman's rank test (GraphPad Prism) was applied to the
154 data to determine if there was a trend in the rank of data, and in some cases a linear regression was used.
155 Discussion of these results is provided in the SI.

156 Overall, these results indicate that all of these assays could benefit from decreasing their
157 variability and increasing the within-laboratory reproducibility. Therefore, the assays were evaluated in
158 greater detail in subsequent sections to discuss possible options for improving their reproducibility and
159 for enabling their use with ENMs.

160 **Evaluation of sources of variability in the assays**

161 Only two steps are similar among all of the assays: the initial cell seeding step and the step related
162 to treating the cells with the ENM and positive chemical controls (see Figure 1), although the ENMs are
163 dispersed in different media among the assays prior to cell treatment. To assess the ENM dispersion, it is
164 important to characterize, preferably using orthogonal techniques, the suspended ENMs such as the
165 extent to which they have agglomerated and their interactions with serum proteins (if in the media). It is
166 important to note though that these measurements are challenging and the results may vary among
167 analytical techniques. Overall, improvements in the repeatability for the ENM dispersion, cell seeding, and
168 dosing steps could help improve the precision of all of the assays. However, the extent that the variability
169 would be decreased would depend on the relative contribution of these steps to the total variability for
170 each assay. It is also evident when viewing the flow charts that there are substantial differences in the
171 number of steps in the assays with the Comet assay and the ELISA method containing the highest number.
172 The instrument used to obtain the final assay readout also varies with two of the assays requiring
173 absorbance measurements, while the other assays utilize flow cytometry, fluorescence or microscopic
174 analysis.

175 The C&E diagrams showed similar sources of variability among assays and key differences (Figures
176 2 through 6). The aspects of the branch colored in orange are designed to highlight differences compared
177 to the previously published C&E diagram for using the MTS assay with ENMs (Figure 2). For all of the assays
178 except for the MTS assay, the positive chemical control to be used to elicit a similar response as the
179 potential impact of the ENM was a significant source of variability. Therefore, this topic will be discussed
180 in depth in subsequent sections.

181 These C&E diagrams were created based upon the specific protocols available in the SI. While
182 there are factors such as the media used or the serum percentage that vary among protocols, these factors
183 are not discussed in depth here since these experimental details are specified in the protocols. Similarly,
184 the use of high throughput testing may help minimize variability in assay results comparing different ENMs
185 if the suspensions are tested at the same time,⁵³ but that is not possible using the protocols described
186 here. In addition, the branches for the C&E diagrams were determined by expert judgement and can be
187 refined based upon robustness testing of the assays⁵⁴ since some of the expected sources of variability
188 may not have a substantial contribution. These sources of variability focus on unintended deviations from
189 the assay protocol (e.g., random variability from pipetting cells among multiple pipetting ejections, or
190 misinterpreting a step of the protocol such as cell rinsing during the MTS assay⁴⁶). The C&E diagrams do
191 not cover mistakes performing the assay protocol (e.g., pipetting double the intended volume or if the
192 cells are contaminated such as with mycoplasma), which should be avoided through a laboratory's quality
193 control system. In other words, to the best of his or her knowledge, the operator executes the protocol
194 correctly, but nevertheless, unnoticed or unavoidable sources of variability lead to varying assay results.

195 To address the sources of uncertainty revealed during the C&E analysis, sensitivity testing can be
196 conducted to quantify the amount of variability from these different sources and control experiments can
197 be performed to evaluate potential artifacts from ENM testing. Overall, potential artifacts from ENM
198 testing depend upon the individual steps in an assay and also the analytical instrument (e.g., plate reader)
199 used to perform the measurement. Thus, there are not overarching control experiments that would be
200 relevant for all four assays. Nevertheless, the conceptual process for identifying control experiments
201 related to artifacts from ENM exposure for these four assays is similar with three primary questions: 1) Is
202 the ENM able to process or modify probe molecule in the absence of cells? 2) Is the ENM fluorescent or
203 absorbent by itself? 3) Is the ENM reducing or increasing an existing signal? These considerations and
204 insights from the C&E diagrams were used to delineate control measurements for each assay several of
205 which are shown in Tables 2 through 5 and are discussed for each assay in the following sections.

206 It is valuable to note that there are key control measurements for instrument performance
207 specific to the instruments used for each of the assays. For example, it is important to ensure that the
208 plate reader is calibrated and to confirm that it provides a linear response for the conditions being used
209 in the assay (e.g., for the fluorescent DCF molecule for the DCF assay). For older plate readers, there is
210 typically only a limited number of settings available, and thus the specific suggested excitation and
211 emission wavelengths for a dye may not be available. Therefore, it is often necessary to determine the
212 optimal settings for each system. It should be noted though that differences in the excitation/emission
213 settings for plate readers can be a source of increased variability among laboratories. There are similar
214 considerations that are relevant for a fluorescent microscope (e.g., linearity, benchmarking) and the flow
215 cytometer (e.g., potential for spillover for the dyes used and assessment of the size of the cells to be
216 tested by measuring their side and forward scattering to determine the binning). Since these aspects of
217 instrument maintenance and calibration are not specific to these assays, they will not be discussed in
218 additional detail in the subsequent sections.

219 In the following sections, each assay will be separately discussed in depth including a general
220 discussion of results from the C&E analysis, control measurements identified from the C&E analysis for
221 the assay (including ENM-specific measurements), and strategies to improve the quality of each assay. It
222 is also important to note that some of these control measurements can be conducted in preliminary
223 experiments such as to identify optimal instrument settings to refine the protocol or to assess for the
224 potential of ENM-induced artifacts. If such artifacts are observed and the impact is substantial, this may
225 indicate that modifications to the protocol are needed or that the assay cannot be used with this ENM. In
226 other cases, modest ENM-induced artifacts may be corrected for using process control measurements,
227 measurements that are made each time the assay is conducted such as the performance of the positive
228 chemical control to provide evidence that the assay is functioning as expected. For example, it may be
229 possible to correct for bias from ENM settling during the MTS assay by performing background subtraction
230 using wells dosed with the same ENM concentration but without cells which undergo all subsequent assay
231 steps.⁴⁶

232 **Detection of reactive oxygen species by H₂DCF-DA assay (DCF assay)**

233 For the DCF assay, a main source of variability relates to the chemical reactivity of the
234 independent positive chemical control (Figure 3). There are no known positive chemical controls for
235 generating consistent quantities of ROS species in this assay, thus hindering the comparability of results
236 with the DCF assay among laboratories or across time in a single laboratory.¹⁴ What is needed is a control
237 substance that can produce a consistent amount of ROS within the cells, is readily available, not cytotoxic,
238 and can be quantified in solution.¹⁴ A positive chemical control that fulfills these criteria could be used as
239 a calibrator in a concentration series to allow the comparability of results across time and space by
240 comparing a response from this positive chemical control such as an EC₅₀ value. Given the inability to find
241 a positive chemical control that can yield a consistent amount of ROS in this assay, the results from this
242 assay are not quantitative and cannot be readily compared across experiments or laboratories. Other
243 unique sources of variability in this assay relate to the instability of the H₂DCF-DA reagent with time and
244 as a result of its potential to be degraded by laboratory light (unpublished results). This can result in
245 variability in the assay results that are challenging to quantify since a positive chemical control capable of
246 generating a consistent quantity of ROS species is not available.

247 There are several key control measurements for the DCF assay in addition to positive chemical
248 control measurements and assessments of potential biases or artifacts from testing ENMs. These control
249 measurements can be assessed during preliminary experiments to assess the robustness of the assay. It
250 is valuable to assess whether loading the dye in the light or dark would influence the assay results,
251 although the H₂DCF-DA molecule, which is first added to the assay, is more stable than the highly light-
252 sensitive fluorescent DCF molecule, which is measured by the assay. In addition, it is also possible to
253 compare the results among H₂DCF-DA provided by different manufacturers to assess to what extent the
254 source of this assay reagent impacts the results. Another consideration is the potential impact of repeated
255 freeze/thaw cycles on the H₂DCF-DA molecules. For each of these topics, it is necessary to keep other
256 factors (e.g., cell handling) as constant as possible or to use cell-free control measurements to assess
257 changes in the H₂DCF-DA molecule. For example, if repeated freeze-thaw cycles are to be tested, one
258 needs to make sure that apparent differences in the H₂DCF-DA with time are not due to changes in the
259 chemical (e.g., Sin-1) that was used to induce ROS which may itself change during storage. As a result of
260 the high Sin-1 reactivity, it is possible that the pipetting time needed to aliquot Sin-1 throughout a whole
261 96-well plate could be sufficiently long to result in changes from the first to the last wells to which Sin-1
262 was added. This can hinder determining if changes in the signal measured are from the Sin-1 or changes
263 to the H₂DCF-DA molecule.

264 There are also several key ENM-specific control experiments for the DCF assay that should be
265 evaluated during preliminary experiments. It is valuable to assess the time needed for ENMs to contact
266 the cells. While the H₂DCF-DA molecule is first allowed to enter the cell after which point the cells cleave
267 the DA and then the excess dye is washed away prior to ENM addition, the time points after the ENM
268 addition should allow sufficient time for the ENM to have contacted the cell before measuring the assay
269 output. It is possible though that the ENMs can produce ROS and impact the deacetylated H₂DCF in the
270 absence of cells, thereby causing a false positive signal. This can be assessed by performing an acellular
271 experiment with the dispersed ENMs and the deacetylated H₂DCF molecule. A positive signal in this
272 acellular control experiment indicates that the ENM can produce ROS extracellularly or that the reactive

273 surface of the ENM can process the dye although not necessarily through the production of ROS. If a
274 positive result in this acellular assay is obtained, a positive response in the cellular assay should be
275 interpreted with caution and additional measurements such as antioxidant biomarkers (e.g., vitamin C or
276 N-acetyl-cysteine) should also be performed before concluding that the ENM causes intracellular ROS. It
277 is also important to assess the extent to which ENMs in the absence of cells and assay reagents can
278 produce results similar to the DCF signal (namely excitation at a wavelength of 485 nm and emission at a
279 wavelength of 528 nm) or can quench the existing fluorescent signal from the DCF molecule^{55,56}. These
280 potential biases can be tested by measuring the fluorescent signal for the dispersed ENMs by themselves
281 or by incubating the fluorescent DCF molecule with dispersed ENMs at a range of ENM concentrations
282 and assessing if there is a change in the fluorescent signal.

283 Numerous steps can be undertaken to further improve the quality of the DCF assay results such
284 as comparing results obtained using reagents from different manufacturers, assessing the impact of
285 freeze-thaw cycles on the H₂DCF-DA, and conducting the ENM related control experiments described
286 above. These ENM related control experiments are critical for proving that the results have not been
287 biased. However, there are limitations regarding the extent to which this assay can be improved given its
288 challenges related to the instability of positive chemical controls which induce ROS thus hindering the
289 ability of results from this assay for different experiments to be compared.¹⁴ Given the sources of
290 uncertainty in this assay, it will be challenging to comprehensively determine that an ENM does not
291 produce ROS above a threshold since the results are currently not comparable among measurements. For
292 example, if ROS is not detected and it is determined that this is not an artifactual result, this finding could
293 stem from multiple factors: insufficient ENM contact with the cells at the time point that ROS is measured
294 given that an increasing fraction of the delivered ENM dose contacts with the cell at longer time points, a
295 lack of ROS production at the time the assay is conducted if the cells have increased their antioxidant to
296 counteract ROS generated by contact with ENMs, the increase in intracellular ROS being too low to detect,
297 or that the ENM does not have the capacity to produce ROS. The highest possible quality for this assay is
298 for there to be grouping of the ROS generating potential perhaps into no detectable ROS generation, weak
299 ROS generation, and strong ROS generation.

300 In addition, this assay is typically conducted in the absence of serum, because the serum can
301 contain esterases which can cleave DA from H₂DCF-DA leaving the molecule vulnerable to extracellular
302 ROS attack thereby producing a false positive signal.^{14,57} Therefore, the results may be challenging to
303 compare to other assays which do use serum because the serum could change the surface of the ENM by
304 producing a protein corona which can impact the ENM interactions with cells, agglomeration behaviors,
305 and *in vitro* dosimetry.^{46,58,59} Nevertheless, it may be possible though to rank the ROS generation of ENMs
306 within a single laboratory if the ENMs are tested during the same experiment or if the results are shown
307 to be comparable across multiple experiments. This assay can also be conducted in two distinct
308 approaches: load the cells with H₂DCF-DA and then apply the ENMs or apply the ENM and then later add
309 the dye. These two approaches can be used to analyze different types of ROS production. If the dye is
310 added first, it is possible to test for an increase in ROS after the first ENM internalization or contact. If the
311 ENM is added first, the initial increase in ROS is not detected, but a longer-term increase in ROS can be
312 measured. However, both of these approaches are very sensitive to kinetics, because the time that the

313 initial increase in ROS occurs depends upon the rate at which the ENMs contact the cells, while longer-
314 term ROS can be mitigated through cellular antioxidant mechanisms. Overall, a positive response from
315 this assay does provide a general indication that ROS are being produced within cells which can be
316 evaluated in more depth if needed with other assays such as more complicated and expensive techniques
317 for acellular measurements (e.g., electron paramagnetic resonance).

318 **Detection of the production and excretion of the pro-inflammatory cytokine IL-8 by ELISA**

319 The ELISA method has several distinct differences from the other assays (Figure 4). For example,
320 it is challenging to assess the purity of the TNF- α positive chemical control. This recombinant protein
321 initiates a cell signaling cascade that results in the cellular production and excretion of IL-8. However, it is
322 not straightforward to assess the purity of TNF- α , and the reagent quality may differ among
323 manufacturers (additional discussion of the positive chemical control is provided in the SI). In addition,
324 other factors such as the quality of the high binding plates can also influence assay results. Similarly, there
325 is substantial variability in the quality of the antibodies used in the assay, but there is not a single
326 universally accepted method for assessing antibody quality.⁶⁰ Antibodies are complex biomolecules that
327 are challenging to fully characterize. Thus, differences among antibodies produced by different
328 manufacturers could influence the assay results. Given that the ELISA protocol requires 22 washing steps,
329 the variability caused by each washing step has the potential to be additive and substantially increase the
330 variability of the assay results. The impact of manual versus automated washing on the variability of assay
331 results can be evaluated during preliminary experiments.

332 One key consideration for this assay is whether to use a commercial kit directly as specified or to
333 titrate antibodies to try to achieve a better sensitivity. For both approaches, it is important to perform a
334 dose-response curve, often known as a standard curve, to evaluate the signal obtained in wells with a
335 known quantity of the cytokine added. The standard curve is needed to assess the performance of the
336 ELISA procedure during that specific experiment (e.g., is the response linear and is the sensitivity
337 sufficient) and to quantify the protein mass in the supernatant of the control and ENM exposed cells. The
338 information from the standard curve can also be used for control charting to evaluate the assay
339 performance across time. It is important to add that the substrate used in the assay (e.g., TMB:
340 tetramethylbenzidine; OPD: o-Phenylenediamine dihydrochloride; or ABTS: Azinobis-
341 ethylbenzothiazoline-6-sulfonic acid-diammonium salt for HRP activation) could result in different
342 sensitivities for the same kit. Several other factors can also impact the assay results including the time
343 allowed for the enzyme-substrate reaction, whether the enzymatic reaction is chemically stopped and the
344 HRP activity (which can change among lots and suppliers and with storage time).

345 Comparable to the DCF assay, it is possible that the ENMs could give a signal similar to the
346 measurand of interest (e.g., absorbance at 630 nm for the measurement of IL-8) or that the ENMs could
347 interact with the assay reagents and decrease the signal strength. This can be evaluated during
348 preliminary experiments for each ENM to be tested. It is possible to conduct acellular control experiments
349 to assess the extent to which ENMs at the cell media concentrations used in the assay provide an
350 absorbance signal at 630 nm. In addition, it is possible to expose the cells to the positive chemical control
351 (typically LPS) for a set exposure duration, remove the supernatant, add ENMs to the supernatant, and

352 then perform the ELISA procedure to assess if the presence of ENMs could cause a bias in the ELISA
353 procedure. Unlike for the DCF assay, another potential bias is for the ENMs to adsorb the IL-8 produced
354 by the cells during the course of the cell exposure, thereby decreasing the measured IL-8 concentration.⁵⁵
355 This can be evaluated by conducting adsorption experiments in the cell media with ENMs at the highest
356 concentration of ENMs to be tested and with a known concentration or range of concentrations for the
357 IL-8 cytokine and then performing the ELISA assay to assess the IL-8 recovery.

358 For ELISA, there have been varying degrees of agreement among the laboratories in the two
359 interlaboratory comparison studies conducted with ENMs.^{24,61} One study showed that all seven
360 laboratories observed significantly increased cytokine levels after cellular exposure to the highest
361 concentration of TiO₂ ENMs as compared to the control samples but the cytokine concentration after TiO₂
362 ENM exposure varied by roughly an order of magnitude.²⁴ In addition, these laboratories showed
363 inconsistent findings about the capacity for certain types of multiwall carbon nanotubes to increase the
364 cytokine concentration. In a second interlaboratory comparison, there were mixed results with two
365 laboratories showing a significantly increased cytokine concentration and two other laboratories showing
366 no effect for cells exposed to silver ENMs.⁶¹ It was unclear why the absolute cytokine concentration varied
367 substantially among laboratories or why there was not better interlaboratory agreement for the second
368 study. However, this variability among laboratories in the absolute cytokine concentration is similar to the
369 intralaboratory variability shown in Figure S2. These findings suggest that the ELISA can potentially yield
370 qualitative interlaboratory agreement (e.g., “yes” or “no” for an increased cytokine concentration
371 compared to the positive chemical control), but it is unclear if this assay can be refined to yield a
372 quantitative agreement. Additional experiments to more carefully evaluate the impact of different
373 sources of variability could be used to further refine the assay and reveal process control measurements
374 to include in the ELISA method to quantify sources of variability more thoroughly during each experiment.
375 However, it is currently unclear which combination of process control measurements would help improve
376 the comparability of results across experiments within a laboratory and among laboratories. It is possible
377 that improving the experimental protocol to decrease the variability in the negative control data could
378 lead to improved comparability, but given the lack of a correlation between the negative and positive
379 control values, it is unclear to what degree this would improve the reproducibility of the assay. Evaluating
380 the robustness of the dose-response relationship for the positive chemical control to changes in the assay
381 protocol may be more valuable.

382 **Detection of cell death mechanism by an Annexin V/PI flow cytometry assay**

383

384 Unlike the ELISA and DCF assay, there are substantially more sources of significant uncertainty in
385 the “instrument performance” branch, largely as a result of challenges related to achieving interlaboratory
386 comparability among flow cytometry results. These challenges are not unique to ENMs. It is challenging
387 with flow cytometry to develop benchmarks to enable comparability among instruments for multiple
388 reasons: it is challenging to develop an absolute calibration; it is challenging to assess the instrument’s
389 linear range of operation; background debris, which is produced by dead cells as they are degraded, may
390 influence the results especially for assays with ENMs, which may have optical properties similar to the

391 dyes used to stain the cells; and ENMs may have unexpected interactions with the cell debris causing
392 measurements potentially similar to apoptotic or necrotic cells. Multiple components of the assay
393 protocol branch may also contribute to the overall uncertainty in the assay results such as the procedure
394 used for preparation of the cells for the flow cytometry analysis including cell harvesting and resuspending
395 the cell pellet, and the quality of the PI and Annexin V reagents given that the Annexin V is only stable in
396 a special buffer. One key step is to determine the gating strategy to be used for all experiments that will
397 be directly compared. The gating strategy and compensation procedures will be especially important for
398 comparisons of results among laboratories. Determining a standardized strategy for gating and
399 compensation is a topic of ongoing research.^{62,63}

400 ENMs may impact specific fluorophores preferentially resulting in either enhanced (false positive)
401 or quenched (false negative) fluorescence signals. For example, Annexin V can be obtained with different
402 fluorophores such as FITC or PE, and "simply" changing the fluorophore may prevent interference
403 reactions. Similar observations have been made for the DNA intercalating dye PI. Changing the "necrosis
404 specific" dye to 7-aminoactinomycin D (7-AAD) eliminated the interference reaction of SiO₂ particles with
405 PI.⁶⁴ The potential for interference of the ENM on the dyes can be evaluated during preliminary
406 experiments. Furthermore, the overall viability status of the cells used for flow cytometry should be
407 assessed by a second independent method such as the trypan blue exclusion assay. This can be conducted
408 each time the assay is performed if feasible. A massive PI signal in flow cytometry without any trypan blue
409 positive cells under the microscope within the same sample indicates ENM interference in at least one of
410 the assays, thereby necessitating further investigations.

411 One important control experiment for this assay is to evaluate if the cell harvesting procedure
412 may produce debris that could then interact with ENMs and hinder the quantification of the apoptotic, or
413 the late apoptotic and necrotic cells. This could be evaluated by a 0 h control experiment in which the
414 cells are immediately harvested after ENM addition and the assay performed. It is possible that ENM
415 agglomeration, ENM interactions with cell debris, or ENM interactions with stained or unstained cells
416 could impact gating choices and subsequent analyses.⁵⁵

417 There have not been any interlaboratory comparisons conducted with flow cytometry to assess
418 necrotic and apoptotic cells after ENM exposure to our knowledge. Nevertheless, several potential topics
419 to refine this assay have been discussed in this manuscript, all of which could support the successful
420 conduct of an interlaboratory comparison. Some of the topics described would be relevant for a range of
421 flow cytometry measurements (e.g., gating), while most of the control experiments relate to ENM-
422 relevant issues such as interference between the ENM signal and that of the measurand. In addition, it is
423 possible to spike in ENMs to the cells obtained during control experiments prior to flow cytometry analysis
424 to evaluate if the presence of ENMs enhances or quenches the signal for apoptotic or necrotic cells.⁶⁴

425 **Detection of cellular DNA damage by the Comet assay**

426 There are numerous components in the "Assay protocol" branch (Branch 5) of the Comet assay
427 C&E diagram that can contribute to uncertainty in assay results (see Figure 6). While some of these are
428 similar to sources of variability determined for other assays (e.g., the cell harvesting step is shared with

429 the flow cytometry assay), many of them relate to the unique aspects of this assay such as the alkaline
430 cell lysis step in the protocol. In addition, the instrument used to perform the gel electrophoresis differs
431 from the instruments used in the other assays. Similar to the other assays, the effectiveness of the positive
432 chemical control is also a major source of variability for the Comet assay.

433 There are numerous control measurements that can be performed to assess the performance of
434 specific aspects of the Comet assay each time the assay is run. There are some commercially available
435 control cells with variable amounts of DNA damage (e.g., leukocytes treated with etoposide⁶⁵ from
436 Trevigen). While these controls are different from positive chemical controls conducted by exposing cells
437 to chemicals during the assay itself, they can provide information about the consistency of the
438 performance of other steps in the assay protocol (e.g., the electrophoresis step). In addition, different
439 dyes to stain DNA after electrophoresis can be utilized. It is important to choose a dye that intercalates
440 into double- and single-stranded DNA and to maintain consistency among the dyes within a study since
441 different dyes may yield varying results.

442 The equipment used to perform the gel electrophoresis can also impact the assay results. For
443 example, it is possible that there is a change in temperature during this step of the assay or that there is
444 heterogeneity in the voltage across the gel. Control measurements could be conducted to evaluate the
445 heterogeneity in the voltage in the gel electrophoresis step by placing slides of control samples at different
446 locations in the tank to see if similar results are achieved. It may also be possible to take steps to control
447 the temperature and pH of the test setup or to measure them before and after the assay is performed to
448 assess if there was a change. These measurements can be performed during preliminary experiments. It
449 may also be helpful to measure the temperature during the gel electrophoresis step each time the
450 experiment is performed to ensure it is consistent among experiments.

451 Another key set of control measurements relates to the microscopic analysis of the comets.⁶⁶ For
452 example, microscopic settings such as the focus and camera exposure time can cause variations in
453 percentage DNA in the tail by up to 40 % although there are steps that can be taken to improve the
454 reproducibility of this component of the assay.⁶⁶ Selection of the comets can also impact the results
455 especially if the slides are not scored blindly (i.e., without knowing which sample the operator is scoring).
456 One approach to minimize this source of variability is to automate the comet selection. It is also possible
457 to analyze the same sample multiple times and potentially using multiple approaches (e.g., manually or
458 with automation, or with using different Comet selection parameters for the automated Comet selection)
459 to assess the variability of this step.

460 There are several relevant potential mechanisms through which ENMs may cause artifacts or
461 biases in the Comet assay: 1) the presence of ENMs may influence the DNA migration rate during the gel
462 electrophoresis step, 2) ENMs may associate with the nucleus, migrate themselves during the gel
463 electrophoresis step, and be misinterpreted as damaged DNA or DNA in the Comet head during the
464 microscopic analysis, and 3) reactive ENMs may damage DNA during the processing steps after the cell
465 exposure (for example through photoactive ENMs being activated by laboratory light during the cell
466 processing steps and causing DNA damage) and then this damage is misinterpreted as having occurred
467 during the cell exposure.^{21,68,69} There are several control measurements that could be performed to

468 evaluate the extent to which these biases or artifacts have occurred. First, a 0-h control experiment could
469 be performed in which ENMs are added to the cells followed by immediate processing of the cells.^{70,71} If
470 there is an increase in apparent DNA damage for the cells with ENMs as compared to untreated cells, this
471 indicates that the ENMs have caused an artifact. It may be possible to wash the nuclei after the treatment
472 process to remove ENMs, but this approach has only recently been tested.⁷² Second, the Comet tail could
473 be analyzed microscopically such as with hyperspectral imaging analysis to evaluate the extent to which
474 ENMs are present in the tail (e.g., using hyperspectral imaging analysis) and the extent to which they could
475 cause apparent DNA damage.⁷³ Third, for potentially photoactive ENMs, it is possible to treat a large
476 number of replicates and then conduct the cell processing for some samples under laboratory light
477 conditions and for other samples with lighting that has a narrow spectrum or reduced light intensity
478 designed to minimize photoactivation.^{21,68,74}

479 One of the main approaches that could be used to improve the quality of the Comet assay is the
480 use of cells with carefully controlled DNA damage such as through irradiation to investigate the impact of
481 different aspects of the assay. The intrinsic variability of biological samples (i.e., cells within a single
482 replicate) typically yields a broad dispersion of results for each replicate. Sample-to-sample variability may
483 also be substantial. Therefore, a batch of frozen, single-use aliquots of consistent cell samples ("reference
484 cells") could, for example, be used to evaluate the impact of different factors with regards to the gel
485 electrophoresis step such as the impact of voltage, temperature, time for unwinding and electrophoresis,
486 and low melting agarose concentration. It could further be used to improve the automation and reliability
487 of the microscopy steps which could produce similar results regardless of the sample operator.

488 A range of potential biases or artifacts have been identified in the Comet assay when used with
489 testing ENMs, but these control experiments and how to handle artifacts when observed can be better
490 elucidated. The extent to which these can be mitigated such as through washing the DNA after the
491 exposure prior to the gel electrophoresis step should be further examined. It is also unclear if it is possible
492 to correct for biases from the ENMs if they are observed or if an alternative genotoxicity assay should be
493 utilized instead. Additional analysis of Comet assay slides with ENMs in the tails could be evaluated as
494 described above to assess the impact of ENMs in the Comet assay. In interlaboratory results from the
495 Comet assay with ENM exposure, there have been mixed results with some laboratories showing non-
496 comparable results.⁷⁵

497 **Conclusions**

498 The steps described in this paper outline a process of evaluating the quality of *in vitro* toxicological
499 assays for use with ENMs using C&E diagrams and describe specific measurements that can be taken to
500 improve the quality of each of the four assays. Producing a C&E diagram is a key step that can be utilized
501 in the development of any robust toxicological assay. The C&E diagram can then guide subsequent
502 robustness testing to quantify different sources of uncertainty. Based on this information, it is possible to
503 design a protocol that includes preliminary control experiments which may need to be performed with
504 each ENM to, for example, identify if potential artifacts are observed, and also process control
505 measurements which yield information about the assay performance each time it is performed. It is
506 important to note that many of the branches in the C&E diagrams are similar among the different assays

507 and were the same as for the MTS assay.¹⁵ Therefore, it may be easier to prepare C&E diagrams for
508 additional assays once it has been conducted for one. In addition, improvements in the precision for
509 shared steps among the assays can result in improvements for all of the assays.

510 The degree to which the precision and robustness of the assay for use with ENMs needs to be
511 improved for a specific application depends upon ensuring that the assay is fit for purpose. For example,
512 assays used for screening purposes may require less precise output and qualitative answers (e.g., yes or
513 no) may be sufficient. For use in quantitative risk assessment or for replacement of *in vivo* assays, the
514 quality of the assay may need to be substantially higher. A promising overall strategy is for some
515 combination of these assays to be used in an integrated approach for testing and assessment (IATA) for
516 use in screening ENMs for additional testing using *in vivo* assays (e.g., for inhalation exposure) or to predict
517 results from *in vivo* assays. It is unlikely that any individual assay will be able to fulfill either of these
518 purposes, but combinations of assays for use in IATAs are more promising. Therefore, it is challenging to
519 predict *a priori* which of these assays will be most helpful for these purposes in the absence of test results
520 and comparisons to a specific set of *in vivo* results.

521 In general, the use of orthogonal methods such as testing cytotoxicity with different approaches
522 is highly valuable for nanocytotoxicity studies and can help identify artifactual results. Nevertheless,
523 increasing the number of assays does also increase the cost and resources required for each ENM to be
524 tested. In addition, it is impossible given the broad range of ENMs that can be synthesized and their
525 variable behaviors and properties to provide descriptive information about which potential biases or
526 artifacts will be the most important for each assay. Previously published results for ENMs with similar
527 properties will likely give guidance about the likelihood of a certain ENM to cause a certain artifact but
528 control experiments need to be conducted on a case-by-case basis. Another key topic for developing
529 assays is their statistical robustness and design to minimize the frequency of false positive and false
530 negative results. Different statistical methods can be utilized for qualitative (e.g., high, low, or no effect)
531 as compared to quantitative results (e.g., an EC₅₀ value with a defined uncertainty) to ensure the assay's
532 statistical robustness.^{76,77} The approach described in this paper can also inform the evaluation of other
533 biological assays for use with ENMs and support their refinement and eventual validation.

534 **Supporting information**

535 Materials and methods section and detailed protocols for all assays, supporting results related to
536 evaluation of the control charting data, supplemental discussion related to positive chemical control
537 considerations, control charting data for all four assays, summary data for different Sin-1
538 concentrations for the DCF assay, and analysis of control charting data for the ELISA method. This
539 information is available free of charge via the Internet at <http://pubs.acs.org>.

540 **Acknowledgements**

541 We acknowledge funding from the NanoScreen Materials Challenge co-funded by the Competence
542 Centre for Materials Science and Technology (CCMX).

543 Certain commercial products or equipment are described in this paper in order to specify adequately the
544 experimental procedure. In no case does such identification imply recommendation or endorsement by
545 the National Institute of Standards and Technology, nor does it imply that it is necessarily the best
546 available for the purpose.

547 **Conflict of interest statement**

548 The authors declare no conflicts of interest.

549 **Corresponding Author:** Elijah J. Petersen, E-mail: elijah.petersen@nist.gov, Phone: 301-975-8142

550

551 **References**

- 552 1. Pavlidou, S.; Papaspyrides, C. D., A review on polymer-layered silicate nanocomposites. *Prog.*
553 *Polymer Sci.* **2008**, *33* (12), 1119-1198.
- 554 2. Petersen, E. J.; Pinto, R. A.; Shi, X. Y.; Huang, Q. G., Impact of size and sorption on degradation of
555 trichloroethylene and polychlorinated biphenyls by nano-scale zerovalent iron. *J. Hazard. Mater.* **2012**,
556 *243*, 73-79.
- 557 3. Petersen, E. J.; Lam, T.; Gorham, J. M.; Scott, K. C.; Long, C. J.; Stanley, D.; Sharma, R.; Liddle, J.
558 A.; Pellegrin, B.; Nguyen, T., Methods to assess the impact of UV irradiation on the surface chemistry
559 and structure of multiwall carbon nanotube epoxy nanocomposites. *Carbon* **2014**, *69*, 194-205.
- 560 4. Potts, J. R.; Dreyer, D. R.; Bielawski, C. W.; Ruoff, R. S., Graphene-based polymer
561 nanocomposites. *Polymer* **2011**, *52* (1), 5-25.
- 562 5. Zou, H.; Wu, S. S.; Shen, J., Polymer/silica nanocomposites: Preparation, characterization,
563 properties, and applications. *Chem. Rev.* **2008**, *108* (9), 3893-3957.
- 564 6. Froggett, S. J.; Clancy, S. F.; Boverhof, D. R.; Canady, R. A., A review and perspective of existing
565 research on the release of nanomaterials from solid nanocomposites. *Part. Fibre Toxicol.* **2014**, *11*.
- 566 7. Kaiser, J. P.; Roesslein, M.; Diener, L.; Wick, P., Human Health Risk of Ingested Nanoparticles
567 That Are Added as Multifunctional Agents to Paints: an In Vitro Study. *Plos One* **2013**, *8* (12).
- 568 8. Kuhlbusch, T. A. J.; Wijnhoven, S. W. P.; Haase, A., Nanomaterial exposures for worker,
569 consumer and the general public. *Nanoimpact* **2018**, *10*, 11-25.

- 570 9. Mitrano, D. M.; Motellier, S.; Clavaguera, S.; Nowack, B., Review of nanomaterial aging and
571 transformations through the life cycle of nano-enhanced products. *Environ. Intl.* **2015**, *77*, 132-147.
- 572 10. Nowack, B.; David, R. M.; Fissan, H.; Morris, H.; Shatkin, J. A.; Stintz, M.; Zepp, R.; Brouwer, D.,
573 Potential release scenarios for carbon nanotubes used in composites. *Environ. Intl.* **2013**, *59*, 1-11.
- 574 11. G. C. Waissi-Leinonen, E. J. Petersen, K. Pakarinen, J. Akkanen, M. T. Leppanen and J. V. K.
575 Kukkonen, Toxicity of fullerene (C60) to sediment-dwelling invertebrate Chironomus riparius larvae,
576 *Environ. Toxicol. Chem.*, 2012, **31**, 2108-2116.
- 577 12. Pakarinen, K.; Petersen, E. J.; Alvila, L.; Waissi-Leinonen, G. C.; Akkanen, J.; Leppänen, M. T.;
578 Kukkonen, J. V. K., A screening study on the fate of fullerenes (nC60) and their toxic implications in
579 natural freshwaters. *Environ. Toxicol. Chem.* **2013**, *32* (6), 1224-1232.
- 580 13. B. C. Nelson, E. J. Petersen, B. J. Marquis, D. H. Atha, J. T. Elliott, D. Cleveland, S. S. Watson, I. H.
581 Tseng, A. Dillon, M. Theodore and J. Jackman, NIST gold nanoparticle reference materials do not induce
582 oxidative DNA damage, *Nanotoxicology*, 2013, **7**, 21-29.
- 583 14. Roesslein, M.; Hirsch, C.; Kaiser, J. P.; Krug, H. F.; Wick, P., Comparability of in Vitro Tests for
584 Bioactive Nanoparticles: A Common Assay to Detect Reactive Oxygen Species as an Example. *Inter. J.*
585 *Mol. Sci.* **2013**, *14* (12), 24320-24337.
- 586 15. Rösslein, M.; Elliott, J. T.; Salit, M.; Petersen, E. J.; Hirsch, C.; Krug, H. F.; Wick, P., Use of Cause-
587 and-Effect Analysis to Design a High-Quality Nanocytotoxicology Assay. *Chem. Res. Toxicol.* **2015**, *28* (1),
588 21-30.
- 589 16. Burden, N.; Aschberger, K.; Chaudhry, Q.; Clift, M. J. D.; Fowler, P.; Johnston, H.; Landsiedel, R.;
590 Rowland, J.; Stone, V.; Doak, S. H., Aligning nanotoxicology with the 3Rs: What is needed to realise the
591 short, medium and long-term opportunities? *Regul. Toxicol. and Pharmacol.* **2017**, *91*, 257-266.
- 592 17. Guggenheim, E. J.; Milani, S.; Röttgermann, P. J. F.; Dusinska, M.; Saout, C.; Salvati, A.; Rädler, J.
593 O.; Lynch, I., Refining in vitro models for nanomaterial exposure to cells and tissues. *NanoImpact* **2018**,
594 *10*, 121-142.
- 595 18. Haase, A.; Lynch, I., Quality in nanosafety - Towards reliable nanomaterial safety assessment.
596 *NanoImpact* **2018**, *11*, 67-68.
- 597 19. Ong, K. J.; MacCormack, T. J.; Clark, R. J.; Ede, J. D.; Ortega, V. A.; Felix, L. C.; Dang, M. K. M.; Ma,
598 G. B.; Fenniri, H.; Veinot, J. G. C.; Goss, G. G., Widespread Nanoparticle-Assay Interference: Implications
599 for Nanotoxicity Testing. *Plos One* **2014**, *9* (3).
- 600 20. Petersen, E. J.; Henry, T. B.; Zhao, J.; MacCusprie, R. I.; Kirschling, T. L.; Dobrovolskaia, M. A.;
601 Hackley, V.; Xing, B.; White, J. C., Identification and Avoidance of Potential Artifacts and
602 Misinterpretations in Nanomaterial Ecotoxicity Measurements. *Environ. Sci. Technol.* **2014**, *48* (8), 4226-
603 4246.
- 604 21. Petersen, E. J.; Reipa, V.; Watson, S. S.; Stanley, D. L.; Rabb, S. A.; Nelson, B. C., DNA Damaging
605 Potential of Photoactivated P25 Titanium Dioxide Nanoparticles. *Chem. Res. Toxicol.* **2014**, *27* (10), 1877-
606 1884.
- 607 22. Riebeling, C.; Piret, J.-P.; Trouiller, B.; Nelissen, I.; Saout, C.; Toussaint, O.; Haase, A., A guide to
608 nanosafety testing: Considerations on cytotoxicity testing in different cell models. *NanoImpact* **2018**, *10*,
609 1-10.
- 610 23. Wang, X.; Duch, M. C.; Mansukhani, N.; Ji, Z.; Liao, Y.-P.; Wang, M.; Zhang, H.; Sun, B.; Chang, C.
611 H.; Li, R.; Lin, S.; Meng, H.; Xia, T.; Hersam, M. C.; Nel, A. E., Use of a Pro-Fibrogenic Mechanism-Based
612 Predictive Toxicological Approach for Tiered Testing and Decision Analysis of Carbonaceous
613 Nanomaterials. *ACS Nano* **2015**, *9* (3), 3032-3043.
- 614 24. Xia, T.; Hamilton, R. F.; Bonner, J. C.; Crandall, E. D.; Elder, A.; Fazlollahi, F.; Girtsman, T. A.; Kim,
615 K.; Mitra, S.; Ntim, S. A.; Orr, G.; Tagmount, M.; Taylor, A. J.; Telesca, D.; Tolic, A.; Vulpe, C. D.; Walker, A.
616 J.; Wang, X.; Witzmann, F. A.; Wu, N.; Xie, Y.; Zink, J. I.; Nel, A.; Holian, A., Interlaboratory Evaluation of

617 in Vitro Cytotoxicity and Inflammatory Responses to Engineered Nanomaterials: The NIEHS Nano GO
618 Consortium. *Environ. Health Perspect.* **2013**, *121* (6), 683-690.

619 25. Collins, A. R.; Annangi, B.; Rubio, L.; Marcos, R.; Dorn, M.; Merker, C.; Estrela-Lopis, I.; Cimpan,
620 M. R.; Ibrahim, M.; Cimpan, E.; Ostermann, M.; Sauter, A.; El Yamani, N.; Shaposhnikov, S.; Chevillard, S.;
621 Paget, V.; Grall, R.; Delic, J.; Goni-de-Cerio, F.; Suarez-Merino, B.; Fessard, V.; Hogeveen, K. N.; Fjellsbo, L.
622 M.; Pran, E. R.; Brzicova, T.; Topinka, J.; Silva, M. J.; Leite, P. E.; Ribeiro, A. R.; Granjeiro, J. M.; Grafstrom,
623 R.; Prina-Mello, A.; Dusinska, M., High throughput toxicity screening and intracellular detection of
624 nanomaterials. *Wiley Interdisc. Rev. Nanomed. Nanobiotechnol.* **2017**, *9* (1).

625 26. Harris, G.; Palosaari, T.; Magdolenova, Z.; Mennecozzi, M.; Gineste, J. M.; Saavedra, L.;
626 Milcamps, A.; Huk, A.; Collins, A. R.; Dusinska, M.; Whelan, M., Iron oxide nanoparticle toxicity testing
627 using high-throughput analysis and high-content imaging. *Nanotoxicology* **2015**, *9*, 87-94.

628 27. Nel, A. E.; Nasser, E.; Godwin, H.; Avery, D.; Bahadori, T.; Bergeson, L.; Beryt, E.; Bonner, J. C.;
629 Boverhof, D.; Carter, J.; Castranova, V.; DeShazo, J. R.; Hussain, S. M.; Kane, A. B.; Klaessig, F.; Kuempel,
630 E.; Lafranconi, M.; Landsiedel, R.; Malloy, T.; Miller, M. B.; Morris, J.; Moss, K.; Oberdorster, G.;
631 Pinkerton, K.; Pleus, R. C.; Shatkin, J. A.; Thomas, R.; Tolaymat, T.; Wang, A.; Wong, J., A Multi-
632 Stakeholder Perspective on the Use of Alternative Test Strategies for Nanomaterial Safety Assessment.
633 *ACS Nano* **2013**, *7* (8), 6422-6433.

634 28. NRC, *Toxicity Testing in the 21st Century: A Vision and a Strategy*. The National Academies Press:
635 2007.

636 29. Godwin, H.; Nameth, C.; Avery, D.; Bergeson, L. L.; Bernard, D.; Beryt, E.; Boyes, W.; Brown, S.;
637 Clippinger, A. J.; Cohen, Y.; Doa, M.; Hendren, C. O.; Holden, P.; Houck, K.; Kane, A. B.; Klaessig, F.;
638 Kodas, T.; Landsiedel, R.; Lynch, I.; Malloy, T.; Miller, M. B.; Muller, J.; Oberdorster, G.; Petersen, E. J.;
639 Pleus, R. C.; Sayre, P.; Stone, V.; Sullivan, K. M.; Tentschert, J.; Wallis, P.; Nel, A. E., Nanomaterial
640 Categorization for Assessing Risk Potential To Facilitate Regulatory Decision-Making. *ACS Nano* **2015**, *9*
641 (4), 3409-3417.

642 30. Clippinger, A. J.; Ahluwalia, A.; Allen, D.; Bonner, J. C.; Casey, W.; Castranova, V.; David, R. M.;
643 Halappanavar, S.; Hotchkiss, J. A.; Jarabek, A. M.; Maier, M.; Polk, W.; Rothen-Rutishauser, B.; Sayes, C.
644 M.; Sayre, P.; Sharma, M.; Stone, V., Expert consensus on an in vitro approach to assess pulmonary
645 fibrogenic potential of aerosolized nanomaterials. *Arch. Toxicol.* **2016**, *90* (7), 1769-1783.

646 31. Landsiedel, R.; Kapp, M. D.; Schulz, M.; Wiench, K.; Oesch, F., Genotoxicity investigations on
647 nanomaterials: Methods, preparation and characterization of test material, potential artifacts and
648 limitations-Many questions, some answers. *Mut. Res.* **2009**, *681* (2-3), 241-258.

649 32. Landsiedel, R.; Sauer, U. G.; Ma-Hock, L.; Schnekenburger, J.; Wiemann, M., Pulmonary toxicity
650 of nanomaterials: a critical comparison of published in vitro assays and in vivo inhalation or instillation
651 studies. *Nanomedicine* **2014**, *9* (16), 2557-2585.

652 33. Krug, H. F.; Wick, P., Nanotoxicology: An Interdisciplinary Challenge. *Ang. Chemi. Intern. Edition*
653 **2011**, *50* (6), 1260-1278.

654 34. OECD. **2018**, Testing Programme of Manufactured Nanomaterials,
655 [http://www.oecd.org/chemicalsafety/nanosafety/testing-programme-manufactured-](http://www.oecd.org/chemicalsafety/nanosafety/testing-programme-manufactured-nanomaterials.htm)
656 [nanomaterials.htm](http://www.oecd.org/chemicalsafety/nanosafety/testing-programme-manufactured-nanomaterials.htm).

657 35. Belyanskaya, L.; Manser, P.; Spohn, P.; Bruinink, A.; Wick, P., The reliability and limits of the MTT
658 reduction assay for carbon nanotubes-cell interaction. *Carbon* **2007**, *45* (13), 2643-2648.

659 36. Guo, L.; Bussche, A. V.; Buechner, M.; Yan, A. H.; Kane, A. B.; Hurt, R. H., Adsorption of essential
660 micronutrients by carbon nanotubes and the implications for nanotoxicity testing. *Small* **2008**, *4* (6),
661 721-727.

662 37. Kroll, A.; Pillukat, M. H.; Hahn, D.; Schnekenburger, J., Interference of engineered nanoparticles
663 with in vitro toxicity assays. *Arch. Toxicol.* **2012**, *86* (7), 1123-1136.

664 38. Worle-Knirsch, J. M.; Pulskamp, K.; Krug, H. F., Oops they did it again! Carbon nanotubes hoax
665 scientists in viability assays. *Nano Lett.* **2006**, *6* (6), 1261-1268.

666 39. Shinohara, N.; Matsumoto, T.; Gamo, M.; Miyauchi, A.; Endo, S.; Yonezawa, Y.; Nakanishi, J., Is
667 Lipid Peroxidation Induced by the Aqueous Suspension of Fullerene C-60 Nanoparticles in the Brains of
668 *Cyprinus carpio*? *Environ. Sci. Technol.* **2009**, *43* (3), 948-953.

669 40. DeLoid, G. M.; Cohen, J. M.; Pyrgiotakis, G.; Demokritou, P., Preparation, characterization, and in
670 vitro dosimetry of dispersed, engineered nanomaterials. *Nat. Prot.* **2017**, *12* (2), 355-371.

671 41. Taurozzi, J. S.; Hackley, V. A.; Wiesner, M. R., Ultrasonic dispersion of nanoparticles for
672 environmental, health and safety assessment - issues and recommendations. *Nanotoxicology* **2011**, *5*
673 (4), 711-729.

674 42. Taurozzi, J. S.; Hackley, V. A.; Wiesner, M. R., A standardised approach for the dispersion of
675 titanium dioxide nanoparticles in biological media. *Nanotoxicology* **2013**, *7* (4), 389-401.

676 43. Drewes, C.; Ojea Jimenez, I.; Mehn, D.; Colpo, P.; Gioria, S.; Bogni, A.; Ponti, J.; Kinsner-
677 Ovaskainen, A.; Gilliland, D.; Riego Sintes, J. *Physicochemical characterisation of gold, silica and silver*
678 *nanoparticles in water and in serum-containing cell culture media, JRC TECHNICAL REPORT*; Publications
679 Office of the European Union, 2017.

680 44. Ellison, S. L. R.; Barwick, V. J., Estimating measurement uncertainty: reconciliation using a cause
681 and effect approach. *Accredit. Qual. Assur.* **1998**, *3* (3), 101-105.

682 45. Food and Drug Administration, Guidance for Industry: Bioanalytical Method Validation Revision
683 1. 2013; p 34.

684 46. Elliott, J. T.; Rosslein, M.; Song, N. W.; Toman, B.; Kinsner-Ovaskainen, A.; Maniratanachote, R.;
685 Salit, M. L.; Petersen, E. J.; Sequeira, F.; Romsos, E. L.; Kim, S. J.; Lee, J.; von Moos, N. R.; Rossi, F.; Hirsch,
686 C.; Krug, H. F.; Suchaoin, W.; Wick, P., Toward Achieving Harmonization in a Nanocytotoxicity Assay
687 Measurement Through an Interlaboratory Comparison Study. *Altex-Alt. Anim. Exper.* **2017**, *34* (2), 201-
688 218.

689 47. Hellack, B.; Nickel, C.; Albrecht, C.; Kuhlbusch, T. A. J.; Boland, S.; Baeza-Squiban, A.; Wohlleben,
690 W.; Schins, R. P. F., Analytical methods to assess the oxidative potential of nanoparticles: a review.
691 *Environ. Sci.: Nano* **2017**, *4* (10), 1920-1934.

692 48. Nel, A.; Xia, T.; Madler, L.; Li, N., Toxic potential of materials at the nanolevel. *Science* **2006**, *311*
693 (5761), 622-627.

694 49. Arts, J. H. E.; Hadi, M.; Irfan, M. A.; Keene, A. M.; Kreiling, R.; Lyon, D.; Maier, M.; Michel, K.;
695 Petry, T.; Sauer, U. G.; Wahrheit, D.; Wiench, K.; Wohlleben, W.; Landsiedel, R., A decision-making
696 framework for the grouping and testing of nanomaterials (DF4nanoGrouping). *Regul. Toxicol.*
697 *Pharmacol.* **2015**, *71* (2), S1-S27.

698 50. Arts, J. H. E.; Irfan, M. A.; Keene, A. M.; Kreiling, R.; Lyon, D.; Maier, M.; Michel, K.; Neubauer, N.;
699 Petry, T.; Sauer, U. G.; Wahrheit, D.; Wiench, K.; Wohleben, W.; Landsiedel, R., Case studies putting the
700 decision-making framework for the grouping and testing of nanomaterials (DF4nanoGrouping) into
701 practice. *Regul. Toxicol. Pharmacol.* **2016**, *76*, 234-261.

702 51. Ayres, J. G.; Borm, P.; Cassee, F. R.; Castranova, V.; Donaldson, K.; Ghio, A.; Harrison, R. M.;
703 Hider, R.; Kelly, F.; Kooter, I. M.; Marano, F.; Maynard, R. L.; Mudway, I.; Nel, A.; Sioutas, C.; Smith, S.;
704 Baeza-Squiban, A.; Cho, A.; Duggan, S.; Froines, J., Evaluating the toxicity of airborne particulate matter
705 and nanoparticles by measuring oxidative stress potential - A workshop report and consensus
706 statement. *Inhal. Toxicol.* **2008**, *20* (1), 75-99.

707 52. Rushton, E. K.; Jiang, J.; Leonard, S. S.; Eberly, S.; Castranova, V.; Biswas, P.; Elder, A.; Han, X. L.;
708 Gelein, R.; Finkelstein, J.; Oberdorster, G., Concept of assessing nanoparticle hazards considering
709 nanoparticle dosimetric and chemical/biological response metrics. *J. Toxicol. Environ. Health-Part A-*
710 *Current Issues* **2010**, *73* (5-6), 445-461.

711 53. El Yamani, N.; Collins, A. R.; Runden-Pran, E.; Fjellsbo, L. M.; Shaposhnikov, S.; Zienolddiny, S.;
712 Dusinska, M., In vitro genotoxicity testing of four reference metal nanomaterials, titanium dioxide, zinc
713 oxide, cerium oxide and silver: towards reliable hazard assessment. *Mutagenesis* **2017**, *32* (1), 117-126.

714 54. Hanna, S. K.; Cooksey, G. A.; Dong, S.; Nelson, B. C.; Mao, L.; Elliott, J. T.; Petersen, E. J.,
715 Feasibility of using a standardized *Caenorhabditis elegans* toxicity test to assess nanomaterial toxicity.
716 *Environ. Sci.: Nano* **2016**, *3* (5), 1080-1089.

717 55. Guadagnini, R.; Kenzaoui, B. H.; Walker, L.; Pojana, G.; Magdolenova, Z.; Bilanicova, D.;
718 Saunders, M.; Juillerat-Jeanneret, L.; Marcomini, A.; Huk, A.; Dusinska, M.; Fjellsbo, L. M.; Marano, F.;
719 Boland, S., Toxicity screenings of nanomaterials: challenges due to interference with assay processes
720 and components of classic in vitro tests. *Nanotoxicology* **2015**, *9*, 13-24.

721 56. Horst, A. M.; Vukanti, R.; Priester, J. H.; Holden, P. A., An Assessment of Fluorescence- and
722 Absorbance-Based Assays to Study Metal-Oxide Nanoparticle ROS Production and Effects on Bacterial
723 Membranes. *Small* **2013**, *9* (9-10), 1753-1764.

724 57. Tarpey, M. M.; Wink, D. A.; Grisham, M. B., Methods for detection of reactive metabolites of
725 oxygen and nitrogen: in vitro and in vivo considerations. *Amer. J. Physiol.-Regulat. Integrat. Compar.*
726 *Physiol.* **2004**, *286* (3), R431-R444.

727 58. Cedervall, T.; Lynch, I.; Lindman, S.; Berggård, T.; Thulin, E.; Nilsson, H.; Dawson, K. A.; Linse, S.,
728 Understanding the nanoparticle–protein corona using methods to quantify exchange rates and affinities
729 of proteins for nanoparticles. *Proceedings of the National Academy of Sciences* **2007**, *104* (7), 2050-
730 2055.

731 59. Lundqvist, M.; Stigler, J.; Elia, G.; Lynch, I.; Cedervall, T.; Dawson, K. A., Nanoparticle size and
732 surface properties determine the protein corona with possible implications for biological impacts. *Proc.*
733 *Nat. Acad. Sci.* **2008**, *105* (38), 14265-14270.

734 60. Bordeaux, J.; Welsh, A. W.; Agarwal, S.; Killiam, E.; Baquero, M. T.; Hanna, J. A.; Anagnostou, V.
735 K.; Rimm, D. L., Antibody validation. *BioTechniques* **2010**, *48* (3), 197-209.

736 61. Piret, J. P.; Bondarenko, O. M.; Boyles, M. S. P.; Himly, M.; Ribeiro, A. R.; Benetti, F.; Smal, C.;
737 Lima, B.; Potthoff, A.; Simion, M.; Dumortier, E.; Leite, P. E. C.; Balottin, L. B.; Granjeiro, J. M.; Ivask, A.;
738 Kahru, A.; Radauer-Preiml, I.; Tischler, U.; Duschl, A.; Saout, C.; Anguissola, S.; Haase, A.; Jacobs, A.;
739 Nelissen, I.; Misra, S. K.; Toussaint, O., Pan-European inter-laboratory studies on a panel of in vitro
740 cytotoxicity and pro-inflammation assays for nanoparticles. *Arch. Toxicol.* **2017**, *91* (6), 2315-2330.

741 62. Glier, H.; Heijnen, I.; Hauwel, M.; Dirks, J.; Quarroz, S.; Lehmann, T.; Rovo, A.; Arn, K.; Matthes,
742 T.; Hogan, C.; Keller, P.; Dudkiewicz, E.; Stüssi, G.; Fernandez, P., Standardization of 8-color flow
743 cytometry across different flow cytometer instruments: A feasibility study in clinical laboratories in
744 Switzerland. *J. Immunol. Meth.* **2017**.

745 63. Wu, D. Y.; Patti-Diaz, L.; Hill, C. G., Development and validation of flow cytometry methods for
746 pharmacodynamic clinical biomarkers. *Bioanalysis* **2010**, *2* (9), 1617-1626.

747 64. Bohmer, N.; Rippl, A.; May, S.; Walter, A.; Heo, M. B.; Kwak, M.; Roesslein, M.; Song, N. W.;
748 Wick, P.; Hirsch, C., Interference of Engineered Nanomaterials in Flow Cytometry: A Case Study. *Coll.*
749 *Surf. B: Biointerfaces* **2018**.

750 65. Hande, K. R., Etoposide: Four decades of development of a topoisomerase II inhibitor. *Eur. J.*
751 *Cancer* **1998**, *34* (10), 1514-1521.

752 66. Braafladt, S.; Reipa, V.; Atha, D. H., The Comet Assay: Automated Imaging Methods for Improved
753 Analysis and Reproducibility. *Sci. Rep.* **2016**, *6*.

754 67. DaNa 2.0, Operating Instructions. (July 20, 2017),
755 <https://nanopartikel.info/en/nanoinfo/methods/992-operating-instructions>


756 68. Karlsson, H. L.; Di Bucchianico, S.; Collins, A. R.; Dusinska, M., Can the Comet Assay be Used
757 Reliably to Detect Nanoparticle-Induced Genotoxicity? *Environ. Mol. Mutag.* **2015**, *56* (2), 82-96.

- 758 69. Rajapakse, K.; Drobne, D.; Kastelec, D.; Marinsek-Logar, R., Experimental evidence of false-
759 positive Comet test results due to TiO₂ particle - assay interactions. *Nanotoxicology* **2013**, *7* (5), 1043-
760 1051.
- 761 70. Lin, M. H.; Hsu, T. S.; Yang, P. M.; Tsai, M. Y.; Perng, T. P.; Lin, L. Y., Comparison of organic and
762 inorganic germanium compounds in cellular radiosensitivity and preparation of germanium
763 nanoparticles as a radiosensitizer. *Inter. J. Rad. Biol.* **2009**, *85* (3), 214-226.
- 764 71. Magdolenova, Z.; Lorenzo, Y.; Collins, A.; Dusinska, M., Can standard genotoxicity tests be
765 applied to nanoparticles? *J. Toxicol. Environ. Health-Part A-Current Issues* **2012**, *75* (13-15), 800-806.
- 766 72. Ferraro, D.; Anselmi-Tamburini, U.; Tredici, I. G.; Ricci, V.; Sommi, P., Overestimation of
767 nanoparticles-induced DNA damage determined by the comet assay. *Nanotoxicology* **2016**, *10* (7), 861-
768 870.
- 769 73. George, J. M.; Magogotya, M.; Vetten, M. A.; Buys, A. V.; Gulumian, M., An Investigation of the
770 Genotoxicity and Interference of Gold Nanoparticles in Commonly Used In Vitro Mutagenicity and
771 Genotoxicity Assays. *Toxicol. Sci.* **2017**, *156* (1), 149-166.
- 772 74. Gerloff, K.; Albrecht, C.; Boots, A. W.; Forster, I.; Schins, R. P. F., Cytotoxicity and oxidative DNA
773 damage by nanoparticles in human intestinal Caco-2 cells. *Nanotoxicology* **2009**, *3* (4), 355-364.
- 774 75. Organization for Economic Cooperation and Development,
775 ENV/JM/MONO(2015)17/PART4/ANN2: Dossier on Titanium Dioxide - Part 4 - NM 102 Annex 2. Paris,
776 France, 2015.
- 777 76. Ellison, S. L. R., Uncertainties in qualitative testing and analysis. *Accredit. Qual. Assur.* **2000**, *5*,
778 346-348.
- 779 77. Ellison, S. L. R.; Gregory, S.; Hardcastle, W. A., Quantifying uncertainty in qualitative analysis.
780 *Analyst* **1998**, *123* (5), 1155-1161.
- 781 78. Drasler, B.; Kucki, M.; Delhaes, F.; Buerki-Thurnherr, T.; Vanhecke, D.; Korejwo, D.; Chortarea, S.;
782 Barosova, H.; Hirsch, C.; Petri-Fink, A.; Rothen-Rutishauser, B.; Wick, P., Single exposure to aerosolized
783 graphene oxide and graphene nanoplatelets did not initiate an acute biological response in a 3D human
784 lung model. *Carbon* **2018**, *137*, 125-135.

785

786


787 **Figures**



788


789 **Figure 1:** Flow charts illustrating steps for MTS, DCF, ELISA, flow cytometry, and Comet assays.

790




791

792 Figure 2: Cause-and-effect diagram of MTS assay. Modified and reprinted with permission from the
 793 American Chemical Society.¹⁵




794

795 Figure 3: Cause-and-effect diagram of the DCF assay. Parts of the diagram in orange font indicate differ
 796 from the cause-and-effect diagram for the MTS assay.




797

798 Figure 4: Cause-and-effect diagram of the ELISA for analyzing IL-8. Parts of the diagram in orange font
 799 indicate differ from the cause-and-effect diagram for the MTS assay.



800

801 Figure 5: Cause-and-effect diagram for the flow cytometry for measuring necrotic and apoptotic cells.
 802 Parts of the diagram in orange font indicate differ from the cause-and-effect diagram for the MTS assay.



803

804 Figure 6: Cause-and-effect for analyzing DNA damage of cells exposed to ENMs with the Comet assay.

805 Parts of the diagram in orange font indicate differ from the cause-and-effect diagram for the MTS assay.

806

807 Table 1: Overview for each of the five assays

Assay	Assay purpose	What is measured
MTS assay	Measures metabolic activity as an indirect estimate of the amount of living cells	The amount of metabolically active cells (i.e., living/healthy cells) is estimated by spectrophotometrically quantifying the amount of MTS reagent metabolized
DCF	Measures reactive oxygen species in cellular systems	A cell-permeable, non-fluorescent dye is added to cells, the cells cleave off the diacetate moiety, the dye is reduced by ROS, and the reduced dye is detected using fluorescence
ELISA Method with IL-8	Measures the cellular release of the cytokine Interleukin-8 (IL-8)	After cellular exposure to a potential stimulus (e.g., ENM), the supernatant is removed and added to a plate which has an antibody attached surface; a second biotinylated antibody is then added to bind to the same cytokine but at a different epitope, and lastly the cytokine concentration is measured spectrophotometrically after horseradish peroxidase linked to avidin reacts using tetramethylbenzidine as the substrate reacts with the biotinylated antibody
Annexin V/ Propidium Iodide Flow Cytometry	Measures cell viability and necrotic/apoptotic cells	Dyes are added to distinguish between viable, necrotic, and apoptotic cells and the cells are measured using flow cytometry
Comet assay	Measures DNA damage	After cell lysis and single cell gel electrophoresis the addition of a DNA-intercalating dye and microscopic imaging allow quantifying DNA breaks in each cell; the relative fluorescence intensity in the tail (% DNA in tail), the tail length as well as the tail moment (taking tail intensity as well as tail length into consideration) are accepted metrics for DNA damage quantification

808

809 Table 2: Control preliminary experiments to assess sources of variability and bias for the DCF assay with
 810 ENMs

Source of variability/bias	Step #	Potential control measurements
Positive chemical control	5	Evaluate potential positive chemical controls ideally to find a compound that always induces the same amount of ROS
Time dependency	17	Perform kinetic measurement because only a small fraction of ENMs may have reached the cells at the bottom of the test plate at early time points (e.g., 4 h). However, it is also valuable to obtain data about the amount of ROS produced shortly after ENM exposure with cells to try to measure the initial release of ROS. It is also possible to quantify ENM uptake in separate experiments prior to the DCF assay (although in practice this may be challenging to measure for some ENMs).
False positive results from extracellular production of ROS	5	Cleave the DA from H ₂ DCF-DA in an acellular experiment and then investigate for potential production of DCF signal during incubation with ENMs as a result of extracellular ROS production. If this is observed, there could be a false positive signal in cellular assays if both extracellular deacetylation and ROS production occur. In addition, if such a signal is observed it could further indicate, that the ENM's reactive surface processes the H ₂ DCF dye without production of ROS.

811

812

813 Table 3: Control preliminary experiments to assess sources of variability and bias for the ELISA protocol
 814 to quantify IL-8 with ENMs

Source of variability/bias	Step #	Potential control measurements
Large number of pipetting/washing steps may increase variability	9, 11, 12, 15	Perform experiments where the washing steps are performed manually or if this step is automated to evaluate the improvement in performance of the assay (e.g., decreased variability of negative control cells and cells exposed to the positive chemical control).
ENM interference with the assay through inherent ENM absorbance or adsorption of assay components	5	Spike ENM into different steps of the assay procedure in the presence and absence of the antigen (used for standard curve preparation) according to SOP on DaNa webpage (V.I.G.O. (2014) "NM interference in an enzyme-linked immunosorbent assay (ELISA) v1.0") ⁶⁷ and evaluate the impact. ⁷⁸
The kinetics of the color reaction may require optimization	17	Take data at several points after substrate application (5, 10, 15, or 20 min) to determine the optimal protocol parameters
Instrument Performance	17	Absorbance standards to characterize function of the instrument; measure with an alternative instrument for comparison
Adsorption of cytokines onto ENMs decreases response	5	A control experiment could be performed to assess adsorption of cytokines by the test ENM. This would consist of incubating the ENMs with a known concentration of IL-8 for a certain interval and then conducting the ELISA method to assess the IL-8 recovery.
Batch to batch variability of ELISA reagents	Not applicable	Run ELISA on the same plate with old and new reagents (standard curve only)

815

816

817

818 Table 4: Control preliminary experiments to assess sources of variability and bias for the Annexin PI flow
 819 cytometry assay with ENMs

Source of variability/bias	Step #	Potential control measurements
Identify a concentration and exposure duration of the positive chemical control to produce a reproducible amount of apoptosis and necrosis cells	5	Conduct preliminary experiments to assess several concentrations and time points. It would be helpful also to test a reference material with a reproducible number of apoptotic and necrotic cells for benchmarking assay performance and comparing to test results.
Interference of ENMs with instrument detection	17	Conduct the following pre-experiments; 1) flow cytometric analysis of cell-free samples (i.e. ENM only) to assess ENMs presence and their interaction with staining solutions; 2) flow cytometric analysis of untreated control cells with ENMs spiked in to assess if ENMs alter the detection of unstained/stained cells or interact with cell debris even if the ENMs are not detected themselves; 3) use other approach to detect fluorescence spectrum of ENMs and ENMs with staining reagents (e.g., test with a plate reader or during physicochemical characterization of the ENMs)

820

821

822 Table 5: Control preliminary experiments to assess sources of variability and bias for the Comet assay
 823 with ENMs

Source of variability/bias	Step #	Potential control measurements
ENM interference with DNA migration	5	Add ENM after cell lysis and assess if there is a change in the Comet size; analyze the Comet tail to assess if ENMs are present such as using electron microscopy or hyperspectral imaging. Assess if washing or separation procedures can be used to minimize ENM concentration in the Comet tail.
Variable voltage during electrophoresis within the plate could influence DNA migration	11	Rotate the slides of control samples to assess if there are differences.
Artifactual light damage after conclusion of cell experiment from photoactive ENMs	5	Conduct the cell exposure as usual and then conduct the remaining steps of the assay with some samples processed with laboratory light and others processed in the dark. Assess if there is a difference between samples processed in the light or dark.

824

825

826 **Author biographies**

827 Elijah J. Petersen: Elijah completed his PhD at the University of Michigan in Environmental Engineering.
828 Then, he completed postdocs at the University of Joensuu (Finland) on a Fulbright scholarship and then
829 the University of Michigan before joining NIST as a National Research Council postdoctoral fellow. He
830 became a staff scientist at NIST in 2010 and works in the Cell Systems Science group in the Biosystems
831 and Biomaterials division. His research currently focuses on the development of robust, reproducible *in*
832 *vitro* test methods. He is an associate editor for Nanotoxicology and Nanoimpact and on the editorial
833 board of Environmental Pollution and Environmental Toxicology and Chemistry.

834 Cordula Hirsch: Cordula Hirsch is a scientist and project leader at Empa developing *in vitro* methods
835 according to the 3Rs principles (replacement, reduction and refinement of animal experiments). A special
836 focus is the improvement and adaptation of *in vitro* cytotoxicity methods for nanosafety as well as
837 nanomedical purposes and to achieve reliable, reproducible and robust test results. She studied biology
838 at the Universities of Konstanz and Freiburg (Germany) and received her PhD in 2007 on cell signaling in
839 the developing and adult nervous system of mice. In 2008, she joined the Particles-Biology interactions
840 laboratory at Empa as a Postdoc.

841 John T. Elliott: John graduated in 1990 with a bachelor's degree in Physics from the University of
842 Massachusetts. After serving as a technical assistant at the Massachusetts General Hospital, he then
843 received a PhD in Physiology and Biophysics from SUNY at Stony Brook in 1999 after which he joined the
844 National Institute of Standards and Technology as a National Research Council postdoctoral fellow. He
845 became a staff scientist at NIST in 2001 and now is the group leader for the Cell Systems Science group
846 in the Material Measurement Laboratory's Biosystems and Biomaterials division.

847 Harald F. Krug: Harald Krug studied Chemistry and Biology. He received his Ph.D. from the Georg-August
848 University in Göttingen for his work in animal physiology. After a postdoc period at the Helmholtz Centre
849 Munich he took over the Department for Environmental Toxicology at the Karlsruhe Institute of
850 Technology. Since 2007, he teaches at the University Berne. At the Swiss Federal Institute for Materials
851 Science & Technology (Empa) he headed the laboratory for Nanomaterials-Biology Interaction and was
852 member of the board from 2010 to 2014. Since 2017 he is retired and is busy with his own company
853 NanoCASE, which he founded in 2014.

854 Leonie Aengenheister: Dr. Leonie Aengenheister received her PhD from ETH Zurich, Switzerland, in 2018
855 for establishing a novel advanced *in vitro* placental model and revealing new insights on placental
856 uptake and translocation of different nanomaterials. She then worked as a postdoctoral researcher at
857 the University of Texas Medical Branch in Galveston, where she designed nanoparticle-based drug
858 carriers destined for therapies during pregnancy. Currently, she continues her investigations on
859 nanomaterial-placenta interactions and potential associated effects on maternal and fetal health as a
860 postdoctoral researcher at Empa - the Swiss Federal Laboratories for Materials Science and Technology
861 in St. Gallen, Switzerland.

862 Ali Talib Arif: Dr. Arif is working as a postdoctoral researcher in the group Environmental Toxicology and
863 Nanotoxicology, Institute for Infection Prevention and Hospital Epidemiology, University Medical Center
864 Freiburg, Germany. He is also a guest scientist at the Kurdistan Institution for Strategic Studies and
865 Scientific Research, Iraq. He holds both B.S and M.Sc. degrees in Biology from the University of
866 Hohenheim, Germany. Dr. Arif's current research studies are deal with molecular mechanisms of
867 environmental chemicals and fine or ultrafine particles toxicity that play a major role to understanding
868 the Health effects associated with particles exposure, so he has published three papers on these topics.

869 Alessia Bogni: Alessia Bogni received a degree in Molecular Biology in 1998 at the University of Milan
870 (Italy), and a PhD in toxicology and pharmacology from University of Konstanz (Germany) in 2003. She
871 gained experience in the field working in Pharmacia&Upjohn (Milan, Italy) and at St. Jude Children
872 Research Hospital (Memphis, TN, USA). She started working at the European Commission Joint Research
873 Centre (Ispra-Italy) in 2006.

874 Agnieszka Kinsner-Ovaskainen: Agnieszka Kinsner-Ovaskainen obtained her M.D. degree from the Medical
875 University of Warsaw, Poland and the Ph.D. degree in Life Sciences from the University of Konstanz,
876 Germany. She joined the European Commission Joint Research Centre (JRC) as official in 2007. Currently
877 she is working in the Health in Society Unit (JRC), responsible for the Central Registries and the European-
878 level coordination of two networks of population-based registries for the surveillance of congenital
879 anomalies (EUROCAT) and cerebral palsy (SCPE).

880 Sarah May: Sarah May studied biological sciences at the University of Konstanz (Germany) and obtained
881 her bachelors degree in 2011, followed by her masters degree in 2014. During her PhD she investigated
882 the effect of different engineered nanomaterials on DNA damage and repair pathways at Empa, St.
883 Gallen in collaboration with the University of Konstanz. In 2018 she received her PhD from the University
884 of Konstanz (Germany). Currently, she is working as medical affairs manager in the medical device
885 industry in Winterthur, Switzerland.

886 Tobias Walser: Dr. Tobias Walser received his PhD in Environmental Engineering on Life Cycle based Risk
887 Assessment of Nanomaterials. After a postdoctoral stay at the U.S. Environmental Protection Agency on
888 modelling the toxic impact of nanomaterials, he worked for the Swiss Federal Office of Public Health.
889 Responsibilities were regulatory risk assessment of nanomaterials, and expert work in various
890 international gremia. In 2017, Tobias Walser founded the company "Vereala", which offers expert
891 knowledge on emerging technologies with a focus on nanomaterials

892 Peter Wick: Peter Wick heads since 2014 the research laboratory for Particles-Biology Interactions at the
893 Federal Laboratories on Materials Science and Technologies Empa in St. Gallen. His general research
894 interest is to study the interactions of nanomaterials with human tissues *in vitro* and *ex vivo* with the
895 purpose to obtain detailed mechanistic understanding about their uptake, accumulation, transport and
896 effect on different types of cells or entire tissue. He is a member of the advisory board of the Swiss Action
897 Plan on Nanomaterials, member of the EDQM working group for NBCs, Editorial Board Member of
898 Nanotoxicology and associated editor of Journal NanoImpact.

899 Matthias Roesslein: Matthias Rösslein studied chemistry at the university of Basel, where he also
900 received his PhD degree in 1989. Afterward he spent to 2 years as a postdoc at the university of Chicago

901 and 4 years as an assistant professor at the physical-chemical institute of the university Zurich, before
902 joining Empa in 1996 as a leading expert in the 'evaluation of measurement uncertainty and metrology'.
903 In 2006 he was appointed a position as 'Senior Scientist' and changed to the "Particles-Biology
904 Interaction" laboratory focusing on the standardization of *in vitro* assays to elucidate the effect of
905 engineered nanomaterials on different cell types.

906