

# Standard Reference Material<sup>®</sup> 2366a

## Cytomegalovirus DNA (Towne<sub>Δ147</sub> BAC) for DNA Measurements

### CERTIFICATE OF ANALYSIS

**Purpose:** The certified values delivered by this Standard Reference Material (SRM) are intended for use in the value assignment of the number of amplifiable genome copies of cytomegalovirus (CMV) per volume sample. When used to value assign practical CMV calibration materials, SRM 2366a can provide traceability to the International System of Units (SI) for measurement results of CMV viral load in tissues or fluids such as plasma.

**Description:** A unit of SRM 2366a consists of one sterile, 0.5 mL perfluoroalkoxy (PFA) fluoropolymer vial containing approximately 150 μL of extracted DNA solubilized in TE buffer, pH 8.0. TE buffer consists of 10 mmol/L 2-amino-2-(hydroxymethyl)-1,3-propanediol (Tris) and 1 mmol/L ethylenediaminetetraacetic acid (EDTA) tetrasodium salt in deionized water. The CMV DNA is in the form of a bacterial artificial chromosome – known as CMV Towne<sub>Δ147</sub> BAC – that contains the genome of the Towne strain of CMV [1,2].

**Certified Values:** The certified values below were determined as the consensus result from five polymerase chain reaction (PCR) assays using a droplet digital (ddPCR) system [3] and a direct measurement of the mean droplet volume. The certified value and its uncertainties were computed using a Gaussian random effect model [4]. A NIST certified value is a value for which NIST has the highest confidence in its accuracy in that all known or suspected sources of bias have been investigated or taken into account [5].

Certified CMV Copy Number Values for SRM 2366a<sup>(a)</sup>

Analyte	Certified Value, $X$	95 % Probability Uncertainty Interval	Effective Standard Uncertainty, $u(X)$	Effective Coefficient of Variation, $CV=100 \times u(X)/X$	Units
CMV copy number	1,796,000	1,746,000 to 1,849,000	26,000	1.5 %	Copies/μL

<sup>(a)</sup> The true value of the analyte in this material is expected to be within the given probability interval with about a 95 % level of confidence. For purposes of uncertainty propagation, use the stated  $u(X)$  or CV, whichever is most convenient for your purpose; these values are associated with 60 degrees of freedom.

**Additional Information:** Additional information is provided in Appendix A.

**Period of Validity:** The certified values delivered by **SRM 2366a** are valid within the measurement uncertainty specified until **01 August 2030**. The certified values are nullified if the material is stored or used improperly, damaged, contaminated, or otherwise modified.

**Maintenance of Certified Values:** NIST will monitor this SRM over the period of its validity. If substantive technical changes occur that affect the certification, NIST will issue an amended certificate through the NIST SRM website (<https://www.nist.gov/srm>). Before making use of any of the values delivered by this material, users should verify they have the most recent version of this documentation, available through the NIST SRM website (<https://www.nist.gov/srm>).

**Safety:** SRM 2366a IS A VIRAL AND BACTERIAL DNA SOURCE MATERIAL. SRM 2366a is a Biosafety Level 1 material and should be handled according to applicable federal, state, and/or local regulations and according to policies and procedures of recipient's organization.

**Storage:** Vials of SRM 2366a should be stored in the dark between 2 °C to 8 °C. **DO NOT FREEZE.**

**Use:** The vial should be mixed by inversion and centrifuged without opening the vial cap prior to removing sample aliquots for analysis. For the certified values to be valid, materials should be withdrawn immediately after opening the vials and processed without delay. Dilutions of these materials may be made as appropriate, but they must be used immediately. Certified values do not apply to any material remaining in recapped vials. **DO NOT EXPOSE ANY DNA SOLUTION TO DIRECT SUNLIGHT.**

**Traceability:** The certified value is traceable to copies per volume, expressed as number of amplifiable CMV entities per microliters of solution, through the confirmation of identity, the ddPCR count measurements, the validity of the Poisson endpoint transformation of counts to entities for digital PCR endpoint assays [6], and calibrated volume measurements made at NIST. Metrological traceability is to the International System of Units (SI) unit for volume.

## REFERENCES

- [1] Wang, W.; Patterson, C.E.; Yang, S.; Zhu, H.; *Coupling Generation of Cytomegalovirus Deletion Mutants and Amplification of Viral BAC Clones*; J. Virol. Meth., Vol. 121, pp. 137–143 (2004).
- [2] GenBank Accession AY315197.2: Human herpesvirus 5 strain Towne, complete genome available at <https://www.ncbi.nlm.nih.gov/nucleotide/AY315197.2> (accessed May 2025).
- [3] Sassenscheidt, J.; Rohayem, J.; Illmer, T.; Bandt, D.; *Detection of B–Herpesviruses in Allogenic Stem Cell Recipients by Quantitative Real-Time PCR*; J. Virol. Meth., Vol. 138, pp. 40–48 (2006).
- [4] Toman, B.; Possolo, A.; *Laboratory Effects Models for Interlaboratory Comparisons*; Accred. Qual. Assur., Vol. 14, pp. 553–563 (2009).
- [5] Beauchamp, C.R.; Camara, J.E.; Carney, J.; Choquette, S.J.; Cole, K.D.; DeRose, P.C.; Duewer, D.L.; Epstein, M.S.; Kline, M.C.; Lippa, K.A.; Lucon, E.; Molloy, J.; Nelson, M.A.; Phinney, K.W.; Polakoski, M.; Possolo, A.; Sander, L.C.; Schiel, J.E.; Sharpless, K.E.; Toman, B.; Winchester, M.R.; Windover, D.; *Metrological Tools for the Reference Materials and Reference Instruments of the NIST Material Measurement Laboratory*; NIST Special Publication (NIST SP) 260-136, 2021 edition; National Institute of Standards and Technology, Gaithersburg, MD (2021); available at <https://nvlpubs.nist.gov/nistpubs/SpecialPublications/NIST.SP.260-136-2021.pdf> (accessed May 2025).
- [6] Vogelstein, B.; Kinzler, K.W.; *Digital PCR*; Proc. Natl. Acad. Sci., Vol. 96, pp. 9236–9241 (1999).
- [7] Sassenscheidt, J.; Rohayem, J.; Illmer, T.; Bandt, D.; *Detection of beta-Herpesviruses in Allogenic Stem Cell Recipients by Quantitative Real-time PCR*; J. Virol. Methods, Vol. 138 (Issues 1–2), pp. 40–48 (2006).
- [8] Bai, X.; Hosler, G.; Rogers, B.B.; Dawson, D.B.; Scheuermann, R.H.; *Quantitative Polymerase Chain Reaction for Human Herpesvirus Diagnosis and Measurement of Epstein-Barr Virus Burden in Posttransplant Lymphoproliferative Disorder*; Clin. Chem., Vol. 43, pp. 1843–1849 (1997).
- [9] Li, H.; Dummer, J.S.; Estes, W.R.; Meng, S.; Wright, P.F.; Tang, Y.-W.; *Measurement of Human Cytomegalovirus Loads by Quantitative Real-Time PCR for Monitoring Clinical Intervention in Transplant Recipients*; J. Clin. Microbiol. Vol. 41, pp. 187–191, (2003).
- [10] Kalpoe, J.S.; Kroes, A.C.M.; de Jong, M.D.; Schinkel, J.; de Brouwer, C.S.; Beersma, M.F.C.; Claas, E.C.J.; *Validation of Clinical Application of Cytomegalovirus Plasma DNA Load Measurement and Definition of Treatment Criteria by Analysis of Correlation to Antigen Detection*; Clin. Microbiol. Vol. 42, pp. 1498–1504 (2004).
- [11] Boeckh, M.; Huang, M.; Ferrenberg, J.; Stevens-Ayers, T.; Stensland, L.; Nichols, W.G.; Corey, L.; *Optimization of Quantitative Detection of Cytomegalovirus DNA in Plasma by Real-Time PCR*; J. Clin. Microbiol., Vol. 42, pp. 1142–1148 (2004).
- [12] ThermoFisher Scientific TaqMan Gene Expression Assay, Cat # 4331182, ID: Pa03453400\_s1 available at [https://www.thermofisher.com/taqman-gene-expression/product/Pa03453400\\_s1?CID=&ICID=&subtype=](https://www.thermofisher.com/taqman-gene-expression/product/Pa03453400_s1?CID=&ICID=&subtype=) (accessed May 2025).

<b>Certificate Revision History:</b> 20 May 2025 (Change of period of validity; updated format; editorial changes); 03 March 2020 (Change of expiration date; editorial changes); 01 September 2015 (Original certificate date).
--

*Certain commercial equipment, instruments, or materials may be identified in this Certificate of Analysis to adequately specify the experimental procedure. Such identification does not imply recommendation or endorsement by the National Institute of Standards and Technology, nor does it imply that the materials or equipment identified are necessarily the best available for the purpose.*

*Users of this SRM should ensure that the Certificate of Analysis in their possession is current. This can be accomplished by contacting the Office of Reference Materials 100 Bureau Drive, Stop 2300, Gaithersburg, MD 20899-2300; telephone (301) 975-2200; e-mail [srminfo@nist.gov](mailto:srminfo@nist.gov); or the Internet at <https://www.nist.gov/srm>.*

**\* \* \* \* \* End of Certificate of Analysis \* \* \* \* \***

# APPENDIX A

**Sample Preparation:** The CMV Towne<sub>Δ147</sub> BAC was obtained from the New Jersey Medical School (H. Zhu, Newark, NJ) and was propagated at NIST in *Escherichia coli* (*E. coli*) cells. Total DNA (CMV and *E. coli*) was isolated, purified, and re-suspended in TE buffer. After concentration adjustment, the material was allowed to equilibrate at 4 °C in a PFA container until vialing. Just prior to vialing, the material was brought to room temperature inside a laminar flow hood and gently mixed. The solution was transferred to sterile PFA vials with an automatic pipette, which were then capped, labeled, and placed into storage at 4 °C.

**Confirmation of CMV Identity:** Sanger and Next Generation Sequencing (NGS) confirm that the CMV DNA in SRM 2366a is completely concordant with reference sequence AY315197.2.

**PCR Assays:** The five quantitative polymerase chain reaction assays used to characterize CMV copy number are listed in Table A1. The primer and probe sequences for these assays are listed in Table A2. Measurements were made using both ddPCR and chamber digital (cdPCR) systems.

Table A1. PCR Assays Used

Name	Locus	Nucleotide position <sup>(a)</sup> (bp)	Amplicon, (bp)	Assay source
CPI	UL54	78,453 to 78,524	72	[7]
gBA	UL56	82,475 to 82,728	254	[8,9]
IE	MIE	170,832 to 170,958	127	[10]
mie11	MIE	172,527 to 172,604	78	Modified from [11]
UL132	UL132	178,461 to 179,522 <sup>(b)</sup>	62	[12]

<sup>(a)</sup> Relative to the CMV reference genome AY315197.2.

<sup>(b)</sup> Commercial assay, position specified relative to the Merlin strain, CMV reference genome AY446894.2.

Table A2. Primer and probe sequences.

Name	Type	Sequence
CPI	Forward primer <sup>(a)</sup>	GGCCGTTACTGTCTGCAGGA
	Reverse Primer <sup>(a)</sup>	GGCCTCGTAGTGAAAATTAATGGT
	Probe – BHQ	CCGTATTGGTGCGGATCTGTTCAA
gBA	Forward primer <sup>(a)</sup>	TACCCCTATCGCGTGTGTTCC
	Reverse Primer <sup>(a)</sup>	ATAGGAGGCGCCACGTATTC
	Probe – BHQ	TTGCTGCCAGCAGATAAGTGGTG
IE	Forward primer <sup>(a)</sup>	CAAGCGCCTCTGATAACCA
	Reverse Primer <sup>(a)</sup>	ACTAGGAGAGCAGACTCTCAGAGGAT
	Probe – BHQ <sup>(a)</sup>	TGCATGAAGGTCTTTGCCAGTACATTCT
mie11	Forward primer	GCATTGTGCTGTGCCTAAGTC
	Reverse Primer	ATGCCTCCAGCGGCTCATG
	Probe – BHQplus	TGGCCTCCACTGTTAGGA
UL132	20X assay – TaqMan MGB	Proprietary sequences

<sup>(a)</sup> Per literature notation, but direction reversed relative to reference sequence AY315197.2.

**Homogeneity Assessment:** Twenty-six units of SRM 2366a were selected randomly using stratified random sampling and the CMV copy number of their DNA solutions assessed with the CPI assay using the ddPCR measurement. The dominant component of observed variance was within-unit repeatability. No systematic trend with production sequence was observed.

**Short-Term Stability Assessment:** Sixteen units of SRM 2366a were selected randomly using stratified random sampling and stored at  $-20\text{ }^{\circ}\text{C}$ ,  $4\text{ }^{\circ}\text{C}$ ,  $21\text{ }^{\circ}\text{C}$ , and  $37\text{ }^{\circ}\text{C}$  (four vials at each temperature) for a period of 30 weeks. CMV copy numbers were assessed ten times during this period with the CP1 assay using both ddPCR and cdPCR measurement systems. The SRM 2366a CMV copy number is thermally stable at  $4\text{ }^{\circ}\text{C}$ ,  $21\text{ }^{\circ}\text{C}$ , and  $37\text{ }^{\circ}\text{C}$  but is unstable to freeze-thaw cycling.

**Certification Measurements:** Seven units of SRM 2366a were selected randomly using stratified random sampling and the CMV copy number of their DNA solutions assessed with the five PCR assays listed in Tables A1 and A2 using both ddPCR and cdPCR measurement systems. Results were obtained from seven measurement campaigns over a period of eight months. Traceable results are currently available only for the ddPCR measurements. However, applying the manufacturer's nominal chamber volume to the cdPCR measurements yields results in excellent agreement with the traceable ddPCR results. Thus, the cdPCR measurements serve to validate the ddPCR results.

\*\*\*\*\* End of Appendix A \*\*\*\*\*