

National Institute of Standards & Technology

Report of Investigation

Reference Material 8852

Ammonium ZSM-5 Zeolite

This Reference Material (RM) is intended to provide a source of ZSM-5 Zeolite for measurement comparisons. The zeolite was donated to the National Institute of Standards and Technology (NIST) by a zeolite-producing company. A unit consists of approximately 35 g to 40 g of material in a sealed container.

Reference and information values derived from chemical, diffraction and particle size measurements are provided in Tables 1 through 8.

Reference Values: Reference values are a best estimate of the true value provided by NIST where all known or suspected sources of bias have not been fully investigated by NIST [1]. Reference values for major components (major elements, loss on ignition [LOI], loss on fusion [LOF]) are provided in Table 1, Si/Al ratio in Table 2, and trace elements in Table 3.

Information Values: Information values are values that may be of interest and use to the user, but insufficient information is available to assess an uncertainty associated with the value [1]. Information values for enthalpy of formation of constituent oxides and liquid water, and elements are provided in Table 4. Information values for unit cell parameters for hydrated ZSM-5 are provided in Table 5. Information values for particle size distribution are provided in Tables 6 and 8; for refractive index determination in Table 7. In addition, spectra for three types of nuclear magnetic resonance (NMR) analyses are provided in Figures 1 through 3; a plot for particle size distribution is provided in Figure 4; electron microscope images are provided in Figure 5 and information concerning the variation in sample mass with change in ambient humidity is provided in Appendix A. Use of this RM in adsorption measurements is provided in Appendix B.

Expiration of Value Assignment: RM 8852 is valid indefinitely, within the measurement uncertainty specified, provided the RM is handled and stored in accordance with instructions given in this Report of Investigation (see "Instructions for Use"). Periodic validation of this RM is not required. This report is nullified if the RM is damaged, contaminated, or otherwise modified.

Maintenance of Reference Values: NIST will monitor this RM over the period of its validity. If substantive technical changes occur that affect the value assignment before the expiration of this report, NIST will notify the purchaser. Registration (see attached sheet or register online) will facilitate notification.

Overall direction and coordination of the analyses were performed at NIST by S. Turner and R.R. Cavanagh.

Analytical measurements were performed at NIST by D.A. Becker, R.A. Fletcher, J.F. Kelly, A.F. Marlow, M.G. Moreno-Ramirez, J.R. Sieber, E.B. Steel, B.H. Toby, S. Turner, J.R. Verkouteren, T.W. Vetter, E.S. Windsor, and R. Zeisler.

NMR studies were performed by G.J. Kennedy and W.S. Borghard at ExxonMobil Research & Engineering (Annandale, NJ). Heat of enthalpy studies were performed by S. Yang and A. Navrotsky at the University of California (Davis, CA).

Statistical consultation for the experimental design and analysis of data was provided by S.D. Leigh, W.F. Guthrie and K.R. Eberhardt at NIST.

Carlos A. Gonzalez, Chief Chemical Sciences Division

Steven J. Choquette, Director Office of Reference Materials

Gaithersburg, Maryland, 20899 Report Issue Date: 13 November 2020 Report Revision History on Page 10 Support aspects involved in the issuance of this RM were coordinated through the NIST Office of Reference Materials.

The zeolite was prepared and packaged by NIST and by Laboratory Quality Services International (LQSI)⁽¹⁾.

INSTRUCTIONS FOR USE

Chemical and Diffraction Analyses: For comparison to chemical and diffraction values, samples should be equilibrated in a constant relative humidity (RH) of 54 $\% \pm 2$ % for at least 48 h. A hydrator (e.g., glove box) containing a saturated solution of Mg(NO₃)₂·6H₂O with excess salt can attain a constant RH of approximately 54 % when the temperature is 20° C. ASTM Standard Practice E 104-02 gives more details on maintaining constant RH by means of aqueous solutions [2]. Samples should be in relatively thin layers (< 0.5 cm) in the hydrator.

If samples are to be weighed for chemical analyses, the samples should be equilibrated to approximately 54 % RH and the samples then weighed in the same ambient humidity. This can be accomplished by weighing samples in a glove box containing a scale. For preparation of samples for diffraction analysis, samples should be equilibrated to approximately 54 % RH in a glove box and then sealed within the glove box. The RH of the glove box should be monitored with a calibrated hygrometer. Samples should not be exposed to ambient laboratory humidity prior to weighing for chemical analyses or during the collection of diffraction data. An example of the variation of sample mass with exposure to a different humidity is given in Appendix A. Refer to section on "Chemical Homogeneity Testing by XRF" for minimum recommended mass for chemical analysis.

Particle Sizing: For particle sizing experiments, special care must be exercised to obtain subsamples that are representative of the particle size distribution of the entire bottle [3–5]. Spin riffling (whole stream rotary sampling) was used to obtain subsamples for characterization of the particle size distribution of this RM. Other methods for subsampling the zeolite material have not been investigated in this study and may result in altered particle size distributions or different uncertainties. In general, the zeolite is highly agglomerated. Therefore, once a representative subsample is obtained, it is recommended that the sample be deagglomerated prior to particle size analysis. Deagglomeration can be achieved by ultrasonification of a suspension of the sample in filtered water. Tests should be run using increasing ultrasonification times to determine the time required to reach a stable particle size distribution. An ice bath may be necessary to reduce the temperature of the suspension during ultrasonification.

The container should be resealed tightly after sampling.

PREPARATION

Sample Preparation: The zeolite material came to NIST in a single metal barrel with a porcelain coating. Samples were scooped out of the original container and placed into aluminized bags using an alternate shovel method. By this method, a scoop was taken from the container and placed in a bag; the next scoop taken from the container was placed in another bag, etc. The procedure was repeated until the original container was empty. This procedure was done to ensure that each bag contained samples from all parts of the original container. Each bag contained less than 9 kg (20 lbs) of material. Sample spin riffling and bottling was contracted out to LQSI, which is a division of Commercial Testing & Engineering Co. (CT&E).

⁽¹⁾ Certain commercial instruments, materials, or processes are identified in this certificate to adequately specify the experimental procedure. Such identification does not imply recommendation or endorsement by the National Institute of Standards and Technology, nor does it imply that the instruments, materials, or processes identified are necessarily the best available for the purpose.

ANALYSES

Methods Used: Chemical analyses were performed by X-ray fluorescence spectrometry (XRF), gravimetry, instrumental neutron activation analysis (INAA) of short, medium and long-lived isotopes, magic angle spinning (MAS) NMR, and calorimetric methods [6].

Unit cell parameters were determined for hydrated samples using synchrotron X-ray diffraction [6].

Particle size analyses were performed using laser light scattering, laser light extinction, and sedimentation methods. In order to obtain representative samples for the analyses [6], a series of mass reduction steps were done to randomly chosen units of material using a spin riffling apparatus. Bottles of zeolite were split into 16 subsamples using a spinning riffler. One of the 16 subsamples was then randomly selected and further split into 8 subsamples using a smaller spinning riffler. An overall riffle split mass reduction of 128 yielded samples of approximately 0.2 g, which is the subsample size used for particle analysis. Tests were performed to determine the amount of sonication necessary to ensure deagglomeration of particles, and samples were ultrasonified prior to particle size measurements. Images of the particles were obtained by scanning transmission electron microscopy (STEM). Electron microscopy samples were not deagglomerated prior to imaging.

Chemical Homogeneity Testing by XRF: Chemical homogeneity was tested for the zeolite for Na, Al, Si, K, Ca, Fe, Zn, and Zr using XRF analysis. No statistically significant differences among bottles were observed for Na, Al, Si, K, Ca, and Fe. A statistical F-test (5 % significance) applied to the data for Zn indicates heterogeneity between bottles at the level of precision of the test data, 0.6 % relative standard deviation for the estimated 0.5 g of material analyzed by the spectrometer. The F-test (5 % significance) applied to the data for Zr indicates heterogeneity between bottles at the level of precision of the test data, 0.8 % relative standard deviation for the estimated 1.0 g of material viewed by the spectrometer. For the major elements, Al and Si, specimens as small as 10 mg will be representative of the bulk material. For the determination of trace elements, the heterogeneity of Zr and Zn must be taken into account and it is recommended that at least 1 g of material be analyzed to obtain results representative of the bulk material.

REFERENCE VALUES: CHEMISTRY

Major Components:

Component	Mass (⁴	s Fra %) ^{(e)(}		Methods
Aluminum $(Al)^{(a)}$ Silicon $(Si)^{(a)}$ Si ^(b) LOI ^{(a)(c)} LOF ^{(a)(d)}	1.396 41.18 45.19 8.50 8.47	± ± ± ±	0.015 0.34 0.13 0.09 0.38	XRF, INAA XRF, gravimetry gravimetry gravimetry XRF

Table 1. Reference Values for Major Components of ZSM-5 Zeolite

^(a) Value relative to the hydrated sample mass.

^(b) Value relative to sample mass ignited at 1000 °C.

- ^(c) Ignited at 1000 °C.
- ^(d) Fused with Li borates at 1100 °C.
- ^(e) The reference values for Al and Si (relative to hydrated mass) are an unweighted mean of the results from two analytical methods. The associated uncertainties are expanded uncertainties about the mean, with coverage factor k = 2 (95 % confidence), calculated by combining a between-method variance [7] with a pooled, within-method variance following the ISO/JCGM and NIST Guides [7]. Where appropriate, components of uncertainty attributable to fit calibration, balance calibration, blank replication, contamination, mechanical loss, and hydration have been included.
- ^(f) The reference values for Si (relative to the ignited sample mass), LOI and LOF are the mean of results obtained by NIST using one analytical method. The expanded uncertainties, U, for Si (relative to ignited sample mass) and LOI are calculated as $U = ku_c$ where u_c is intended to represent, at the level of one standard deviation, the combined standard uncertainty calculated according to the ISO and NIST Guides [8]. The coverage factor, k, is determined from the Student's t-distribution with the appropriate associated degrees of freedom for 95 % confidence for each analyte. The expanded uncertainty for LOF is calculated as $U = ku_c$ where u_c represents the standard uncertainty of the mean for the NIST values, and k = 2.131 is the Student's t multiplier corresponding to 95 % confidence for 15 determinations.

Table 2. Reference Value for Atomic Ratio of Si (Relative to Hydrated Sample Mass) to Al^(a)

Elements	Ratio	Methods
Si/Al	28.34 ± 0.39	XRF, INAA, gravimetry

^(a) The reference atomic ratio value for Si/Al is calculated using the reference values for the elements from Table 1 and the atomic mass of the elements. The corresponding uncertainties are determined by propagation-of-error [9] for ratios from the unexpanded uncertainties of Table 1, and then expanded by a coverage factor k = 2 (95 % confidence).

Trace Elements:

Element ^(a)	Mass Fraction (mg/kg) ^(b)		
Iron (Fe)	222.6	±	9.4
Cerium (Ce)	1.708	±	0.051
Cobalt (Co)	0.1848	±	0.0051
Chromium (Cr)	4.48	±	0.13
Hafnium (Hf)	1.337	±	0.045
Lanthanum (La)	0.471	±	0.014
Antimony (Sb)	1.581	±	0.045
Scandium (Sc)	0.0943	±	0.0030
Samarium (Sm)	0.0815	±	0.0025
Thorium (Th)	0.335	±	0.0098
Tungsten (W)	2.240	±	0.068
Zinc (Zn)	0.749	±	0.048

- ^(a) Trace element analyses by INAA were conducted for Ag, As, Au, Cd, Ce, Co, Cr, Cs, Fe, Hf, La, Mo, Rb, Sb, Sc, Se, Sm, Sn, Th, W, and Zn. Several elements assayed were at or below the detection limit achieved for this material with the procedure used. These elements and their detection limits are: Ag (15 μg/kg), As (10 μg/kg), Au (0.2 μg/kg), Cd (0.3 mg/kg), Cs (10 μg/kg), Mo (0.6 mg/kg), Rb (0.1 mg/kg), Se (25 μg/kg), and Sn (2.5 mg/kg). Analysis for Na was conducted by XRF. Na was not detected above the detection limit for the method of 40 mg/kg.
- ^(b) The reference value for Fe (relative to hydrated mass) is an unweighted mean of the results from two analytical methods (XRF and INAA). The uncertainty listed with the value is an expanded uncertainty about the mean, with coverage factor k = 2 (95 % confidence), calculated by combining a between-method variance [7] with a pooled, within-method variance following the ISO/JCGM and NIST Guides [8]. The reference values for trace elements other than Fe are the mean of results obtained by NIST using one analytical method (INAA applied to 16 samples). The expanded uncertainties, U, are calculated as $U = ku_c$ where u_c is intended to represent, at the level of one standard deviation, the combined standard uncertainty calculated according to the ISO/JCGM and NIST Guides [8]. The coverage factor, k, is determined from the Student's *t*-distribution with the appropriate associated degrees of freedom for 95 % confidence for each analyte.

Comment 1: The compositional formula of RM 8852 determined from chemical analyses is:

 (NH_4) is presumed to provide charge balance for the Al content. The LOI value in Table 1 is presumed to include both H₂O and NH₃. No crystalline impurities were noted in the diffraction work. The uncertainties for the elemental and water content in the formula for ZSM-5 are: Al \pm 0.04, Si \pm 0.8, H₂O \pm 0.3. The uncertainties are propagated from the uncertainties in mass fractions given in Table 1. Trace elements have not been included in the formula.

INFORMATION VALUES: CHEMISTRY

NMR Analyses for Si Coordination: The ²⁹Si MAS NMR spectra obtained from three samples consisted of two main peaks due to Si(OSi)₄ and Si(OAl)(OSi)₃.




Figure 1. ²⁹Si MAS NMR spectrum obtained from ZSM-5 Zeolite.

NMR Analyses for Al Coordination: The ²⁷Al MAS NMR spectra obtained from three samples exhibited a single resonance at 54.9 ppm corresponding to tetrahedral framework Al. There is no evidence for nonframework, octahedrally coordinated aluminum in this sample.




Figure 2. ²⁷Al MAS NMR spectrum obtained from ZSM-5 Zeolite.

NMR Analyses for H in Proximity to Si: An example ²⁹Si cross polarization (CP) MAS NMR spectrum is shown in Figure 3. The CP experiment enhances the signals from those silicon nuclei that are in close proximity to hydrogen, such as covalently bound silanols or trapped hydrocarbons. The signal to noise ratio of a CP-MAS NMR spectrum when compared to that of the ²⁹Si MAS NMR experiment is a qualitative indicator of the presence of these species.




Figure 3. ²⁹Si CP-MAS NMR spectrum obtained from ZSM-5 Zeolite.

Enthalpy of Formation:

Table 4. Enthalpy of Formation from Constituent Oxides and Liquid Water, and Constituent Elements at 25° C

Constituents	Enthalpy of Formation (kJ/mol) ^(a)
Oxides and Liquid Water Elements	$-6.26 \pm 0.54 \\ -990.8 \pm 1.1$

^(a) The uncertainties are propagated from calorimetric data and from literature values for the constituent oxides or elements (coverage factor k = 2,95 % confidence).

INFORMATION VALUES: DIFFRACTION

Synchrotron X-ray Diffraction of Hydrated Samples:

Table 5. Unit Cell Parameter of Hydrated Ammonium ZSM-5 Zeolite (Space Group Pnma)

Axis	Value (nm) ^(a)		
<u>a</u> <u>b</u> <u>c</u>	$\begin{array}{rrrrrrrrrrrrrrrrrrrrrrrrrrrrrrrrrrrr$		

^(a) The uncertainty is one standard deviation based on 16 samples.

Comment 2: The unit cell parameters are for information only. The uncertainties for the reported data do not include all sources of bias. For this material, a possible trimodal distribution of unit cell values were obtained.

INFORMATION VALUES: PARTICLE SIZING

The distributions reported for the particle size measurement methods reported below should be considered valid only for the specific method used. Variations in values among methods may result from effects of particle characteristics such as morphology.

Laser Light Scattering:

Table 6. Cumulative Volume Percent Less than Particle Diameter of Ten Samples^(a)

Cumulative Volume (%)	Mean Diameter (µm)	Standard Deviation of the Mean (µm)
5	0.231	0.002
10	0.268	0.002
15	0.298	0.001
20	0.324	0.001
25	0.350	0.001
30	0.375	0.001
35	0.400	0.001
40	0.427	0.001
45	0.454	0.001
50	0.484	0.001
55	0.517	0.002
60	0.553	0.002
65	0.595	0.003
70	0.646	0.004
75	0.712	0.005
80	0.803	0.007
85	0.971	0.012
90	1.698	0.043
95	3.787	0.058

^(a) Measured by the laser light scattering method.




Figure 4. Distribution of differential volume fraction (%) vs. particle diameter as measured by the laser light scattering method. Each curve represents the analysis for a different bottle of the five bottles analyzed.

Table 7. Refractive Index Determination

Refractive Index^(a)

1.472 to 1.474

^(a) A value of 1.47 was used in Mie modeling of the particle size distribution for the laser light scattering method since values to 2 decimal places are considered sufficient for this purpose.

Laser Light Extinction: The number of particles observed by this method is low and indicates that the bulk of the particles are less than 1 µm in diameter. This is consistent with results for laser light scattering method.

X-ray Sedimentation:

Cumulative Volume (%)	X-ray Sedimentation Mean Diameter (µm) ^(b)	Laser Light Scattering Mean Diameter (µm) ^(c)
10	0.21	0.268 ± 0.002
25	0.30	0.350 ± 0.001
50	0.45	0.484 ± 0.001
75	0.74	0.712 ± 0.005
90	1.41	1.698 ± 0.043

Table 8. Information Values Obtained by X-ray Sedimentation^(a)

^(a) Values from another method for comparison

^(b) Values for X-ray sedimentation are included for particle diameters less than 0.5 μm but the values are less reliable than those from larger particles due to Brownian motion.

^(c) Uncertainty values represent the standard deviation of the mean.




Figure 5. Dark field scanning transmission electron microscopy images of ZSM-5 Zeolite.

REFERENCES

- [1] May, W.; Parris, R.; Beck, C.; Fassett, J.; Greenberg, R.; Guenther, F.; Kramer, G.; Wise, S.; Gills, T.; Colbert, J.; Gettings, R.; MacDonald, B.; *Definitions of Terms and Modes Used at NIST for Value-Assignment of Reference Materials for Chemical Measurements*; NIST Special Publication 260-136; U.S. Government Printing Office: Washington, DC (2000); available at https://www.nist.gov/system/files/documents/srm/SP260-136.PDF (accessed Nov 2020).
- [2] ASTM E104-02; *Standard Practice for Maintaining Constant Relative Humidity by Means of Aqueous Solutions*; Annu. Book ASTM Stand., Vol. 11.03 (2012).
- [3] Kaye, B.H.; *Investigation into the Possibilities of Developing a Rate Method of Sieve Analysis*; Powder Metall., Vol 5, pp. 213–234 (1962).
- [4] Allen, T.; Khan, A.A.; *Critical Evaluation of Powder Sampling Procedures*; Chem. Engr., Vol. 238, pp. CE108–CE112 (1970).
- [5] British Standard 3406-1; *Methods for Determination of Particle Size Distribution Part 1: Guide to Powder Sampling;* British Standards Institution, London UK (1986).
- [6] Turner, S.; Sieber, J.R.; Vetter, T.W.; Zeisler, R.; Marlow, A.F.; Moreno-Rameriz, M.G.; Davis, M.E.; Kennedy, G.J.; Borghard, W.G.; Yang, S.; Navrotsky, A.; Toby, B.H.; Kelly, J.F.; Fletcher, R.A.; Windsor, E.S.; Verkouteren, J.R.; Leigh, S.D.; *Characterization of Chemical Properties, Unit Cell Parameters and Particle Size Distribution of Three Zeolite Reference Materials: RM 8850 Zeolite Y, RM 8851 Zeolite A and RM 8852 Ammonium ZSM-5 Zeolite;* Micropor. Mesopor. Mater., Vol. 107, pp. 252–267 (2008).
- [7] Levenson, M.S.; Banks, D.L.; Eberhardt, K.R.; Gill, L.M.; Guthrie, W.F.; Liu, H.K.; Vangel, M.G.; Yen, J.H.; Zhang, N.F.; An Approach to Combining Results From Multiple Methods Motivated by the ISO GUM; J. Res. Natl. Inst. Stand. Technol., Vol. 105, pp. 571–579 (2000); available online at https://www.ncbi.nlm.nih.gov/pmc/articles/PMC4877152/ (accessed Nov 2020).
- [8] JCGM 100:2008; Evaluation of Measurement Data Guide to the Expression of Uncertainty in Measurement (GUM 1995 with Minor Corrections); Joint Committee for Guides in Metrology (JCGM) (2008); available at https://www.bipm.org/utils/common/documents/jcgm/JCGM_100_2008_E.pdf (accessed Nov 2020); see also Taylor, B.N.; Kuyatt, C.E.; Guidelines for Evaluating and Expressing the Uncertainty of NIST Measurement Results; NIST Technical Note 1297; U.S. Government Printing Office: Washington, DC (1994); available at https://www.nist.gov/pml/nist-technical-note-1297 (accessed Nov 2020).
- [9] Ku, H.H.; Notes on the Use of Propagation of Error Formulas; J. Res. Natl. Bur. Stand., Vol. 70C, pp. 263–273 (1966); available online at https://nvlpubs.nist.gov/nistpubs/jres/70C/jresv70Cn4p263_A1b.pdf (accessed Nov 2020).
- [10] Nguyen, H.G.T.; Espinal, L.; van Zee, R.D.; Thommes, M.; Toman, B.; Hudson, M.S.L.; Mangano, E.; Brandani, S. Broom, D.P.; Benham, M.J.; Cychosz, K.; Bertier, P.; Yang, F.; Krooss, B.M.; Siegelman;, R.L.; Hakuman, M.; Nakai, K.; Ebner, A.D.; Erden, L.; Ritter, J.A.; Moran, A.; Talu, O.; Huang, Y.; Walton, K.S.; Billemont, P.; De Weireld, G.; *A Reference High-Pressure CO₂ Adsorption Isotherm for Ammonium ZSM-5 Zeolite: Results of an Interlaboratory Study*; Adsorption, Vol. 24, pp. 531–539 (2018).
- [11] Span, R.; Wagner, W; A new equation of state for carbon dioxide covering the fluid region from the triple-point temperature to 1100 K at pressures up to 800 MPa; J. Phys. Chem. Ref. Data, Vol. 25, pp. 1509-1597 (1996).
- [12] NIST/ARPA-E; *Database of Novel and Emerging Adsorbent Materials*, National Institute of Standards and Technology, Gaithersburg, Maryland USA (2014); available at https://adsorption.nist.gov/isodb/index.php#home (accessed Nov 2020).

Report Revision History: 13 November 2020 (Correction to reference numbers is Appendix B; editorial changes); 28 September 2020 (Change of expiration date; addition of Appendix B; clarification for LOF Method in Table 1; editorial changes); 26 July 2006 (Original Report of Investigation).

Users of this RM should ensure that the Report of Investigation in their possession is current. This can be accomplished by contacting the SRM Program: telephone (301) 975-2200; fax (301) 948-3730; e-mail srminfo@nist.gov; or via the Internet at https://www.nist.gov/srm.

APPENDIX A

VARIATION IN SAMPLE MASS WITH CHANGE IN AMBIENT HUMIDITY

Experiments were conducted to determine the magnitude of variation in sample mass with change in ambient humidity. One of the experiments is described and its results given as information values so that users of the RM will be aware of the magnitude of mass change with exposure to two relative humidities ($33 \% \pm 2 \%$ RH and $54 \% \pm 2 \%$ RH).

Approximately 0.85 g of zeolite was placed in an aluminum weighing dish. The sample was placed in a hydrator (in this case, a glass desiccator containing a saturated solution of MgCl₂·6H₂O with excess salt) and allowed to equilibrate at 33 % \pm 2 % RH over several days. The hydrator with the sample was placed in a glove box containing a weighing scale and several pans of MgCl₂·6H₂O salt in solution. After the glove box reached 33 % \pm 2 % RH, the sample was removed from the hydrator and the sample mass was determined. The sample was placed back in the hydrator. The pans of MgCl₂·6H₂O salt solution in the glove box were exchanged for several pans containing Mg(NO₃)₂·6H₂O salt solutions. After the glove box was stable at 54 % \pm 2 % RH, the sample was removed from the hydrator he glove box was stable at 54 % \pm 2 % RH, the sample was removed from the first hour and periodically over the next several days.

The change in mass for approximately the first hour is shown in Figure A1. The results show a mass difference of approximately 1.2 % between the 33 % \pm 2 % RH and 54 % \pm 2 % RH provided by the two salt solutions. The greatest change in mass occurred in the initial few minutes of exposure to 54 % \pm 2 % RH. After 48 h, the mass gain stabilized at approximately 1.25 %.




Figure A1. Variation in mass of ZSM-5 Zeolite as the ambient humidity was changed from 33 $\% \pm 2 \%$ RH (time = 0) to 54 $\% \pm 2 \%$ RH (time > 0 min).

Page 11 of 12

APPENDIX B

USE OF NIST REFERENCE MATERIAL 8852 IN ADSORPTION MEASUREMENTS

Reference Material (RM) 8852 was used in an interlaboratory study to develop a reference, high-pressure carbon dioxide surface excess adsorption isotherm [10]. This Appendix summarizes the procedure for realizing this isotherm. Full details are given in the Supplemental Material to Reference [10].

Key Measurement Steps

- The sample mass to be used should be optimized to achieve the best signal-to-noise ratio, according to the experimental parameters of the instrument being used.
- The adsorptive, carbon dioxide, should have a purity of not less than 99.999 %.
- The zeolite should be outgassed by heating slowly to 350 °C (1 °C per minute, nominally) and held at that temperature for at least 12 hours under high vacuum (<1 cPa) with continuous pumping. If the outgassing is performed in a separate manifold, prevent contact with air when transferring to the analysis instrument.
- Carbon dioxide adsorption isotherms should be measured at 20 °C.
- The same measurement is also to be done without the adsorbent to perform a blank correction. This blank should be subtracted from the adsorption isotherm [10].
- The Span and Wagner equation of state [11] should be used to determine the density/compressibility of the gas for the determination of surface excess uptake.

The isotherm is represented by the empirical function

$$n_{ex} = d / \left[1 + \exp\left(\left(-\ln\left(P\right) + a\right) / b\right) \right]^c \pm U_{k=2},$$

where n_{ex} is surface excess amount-of-substance adsorbed per unit mass (mmol/g), *P* is equilibrium pressure (MPa); a = -6.22, b = 1.97, c = 4.73, and d = 3.87; and $U_{k=2} = 0.075$ mmol/g provides the k = 2 prediction interval. The function is valid from 1 kPa to 4.5 MPa. The reference isotherm is also available in a tabular form on-line [12]. Figure B1 shows the reference isotherm with its associated prediction interval, along with a measured isotherm.




Figure B1: Surface excess adsorption (n_{ex}) as function of pressure (P) for RM 8852.